

HEALTHCARE FACILITIES HAVE AN OPEN 30 DAY ACCOUNT

NO MINIMUM ORDER - PURCHASE ORDERS ACCEPTED

F.O.B. Origin. • Prices subject to change without notice.
 25% restocking charge prior to 30 days. • Please call for a RGA Number.
 Products may NOT be returned after 30 days.

PRICING ON WEBSITE SUPERSEDES CATALOG PRICES

INDEX





ADAPTERS


Disposable 19
 DISS 27
 Female Couplers 25
 Male Quick Connects 24
 NPT Fittings 22
 Y Blocks and T Blocks 23
 BAGS, EQUIPMENT COVERS 29
 COMPRESSORS 30
 CYLINDER CARTS AND STANDS 9 - 11
 FITTINGS (under ADAPTERS)
 FLOWMETERS 5 - 8
 FLOWMETER ACCESSORIES 31


HOSES 17, 18
 MISCELLANEOUS 29, 30
 OXYGEN ANALYZERS AND SENSORS 28
 RAIL SYSTEMS 16
 REGULATORS
 Air 15
 CO2 15
 Heliox 15
 Oxygen 12, 13
 Oxygen - Specialty 14
 Suction 2 - 4
 RESPIROMETERS 21

SUCTION REGULATOR SPECIALS LIMITED QUANTITY

				
Ohio Medical Continuous/Intermittent Suction Regulator with Trap	Precision Medical 2 Mode 0-200 mm Hg Continuous Suction Regulator	Allied Healthcare 2 Mode 0-300 mm Hg Continuous Suction Regulator	Ohio Medical Continuous "Push to Set" 3 Mode Suction Regulator with trap	Allied Healthcare 3 Mode 0-300 mm Hg Continuous Suction Regulator




	Ohmeda Wall Adapter	Chemtron Wall Adapter	Puritan Wall Adapter	DISS Wall Adapter
				
	Part #	Part #	Part #	Part #
Ohio Medical Intermittent	LQ-45	LQ-47	LQ-48	LQ-46
Precision	LQ-13	LQ-14	LQ-15	LQ-16
Allied 2 Mode	LQ-41	LQ-42	LQ-43	LQ-44
Ohio Medical Continuous	LQ-49	LQ-51	LQ-52	LQ-50
Allied 3 Mode	LQ-58	LQ-60	LQ-57	LQ-59

 All suction regulators above come standard with tubing nipple outlet.

 For DISS Male outlet, add a "D" to end of part number.

Toll Free Phone
 888-431-2727
 1-860-871-4156

Toll Free Fax
 888-337-0443
 1-860-870-0443

SUCTION REGULATORS – by Precision Medical
CONTINUOUS AND CONTINUOUS/INTERMITTENT

NEW LOWER PRICING

MADE IN USA
10 YEAR WARRANTY

USER FRIENDLY LARGE KNOB
EASY TO READ FACE

2 Mode Continuous



Wall Adapter
Inlet

Tubing Nipple
Outlet

DISS Male
Outlet

DISS Female
Ohmeda
Chemetron
Puritan



PM3003
PM3004
PM3005
PM3008



PM3014
PM3015
PM3016
PM3019

3 Mode Continuous/Full line



Wall Adapter
Inlet

Tubing Nipple
Outlet

DISS Male
Outlet

DISS Female
Ohmeda
Chemetron
Puritan



PM3103
PM3104
PM3105
PM3108



PM3114
PM3115
PM3116
PM3119

3 Mode Continuous/Intermittent



Wall Adapter
Inlet

Tubing Nipple
Outlet

DISS Male
Outlet

DISS Female
Ohmeda
Chemetron
Puritan



PM3303
PM3304
PM3305
PM3308



PM3314
PM3315
PM3316
PM3319

WALL ADAPTERS (INLET)



DISS Female
w/hand nut



Ohio (Ohmeda)




Chemetron (NCG)




Puritan

VACUTRON SUCTION REGULATORS by Allied

The new, sleekly contoured Vacutron® is a precision engineered instrument that combines the most desirable features for suction regulators. Smaller in size, but big on benefits, it offers advantages that include:



Collection Canisters



Disposable Collection Canisters

Size	Outlet	Units/Case	Part Number	Price
2400ML	DISS	36/CASE	20-08-0001	
2400ML	STEM	36/CASE	20-08-0002	
1500ML	STEM	48/CASE	20-08-0003	
1500ML	DISS	48/CASE	20-08-0004	
2400ML	DISS	36/CASE w/tubing	20-08-0005	
2400ML	STEM	36/CASE w/tubing	20-08-0006	
1500ML	STEM	42/CASE w/tubing	20-08-0007	
1500ML	DISS	42/CASE w/tubing	20-08-0008	



Diss Female w/hand nut



Chemetron (NCG)



Ohio (Ohmeda)



Puritan



Tubing Nipple



DISS Male

Continuous Suction Regulator

Part #	Inlet	Outlet	Price
22-13-1106	Ohmeda	DISS Male	
22-13-1107	DISS Female	DISS Male	
22-13-1108	Chemetron (NCG)	DISS Male	
22-13-1105	Puritan-Bennett	DISS Male	
22-13-1206	Ohmeda	Tubing Nipple	
22-13-1207	DISS Female	Tubing Nipple	
22-13-1208	Chemetron (NCG)	Tubing Nipple	
22-13-1205	Puritan-Bennett	Tubing Nipple	

Continuous/Intermittent Suction Regulator

22-15-1106	Ohmeda	DISS Male	
22-15-1107	DISS Female	DISS Male	
22-15-1108	Chemetron (NCG)	DISS Male	
22-15-1105	Puritan-Bennett	DISS Male	
22-15-1206	Ohmeda	Tubing Nipple	
22-15-1207	DISS Female	Tubing Nipple	
22-15-1208	Chemetron (NCG)	Tubing Nipple	
22-15-1205	Puritan-Bennett	Tubing Nipple	

Suction Canister Holder and Wall Bracket



White Plastic Holder
20-02-0150



Wall Bracket
22-11-0006

For suction accessories please visit our website at:

www.rcmedical.com

Tubing
Fittings
Vac Trap

Stands
Bacteria Filters
Wye Blocks

PRECISION MEDICAL CHROME BODY FLOWMETER

**5 Year
Warranty**



4MFA1001



1MFA9001



3MFA1001



0 - 1 LPM PEDIATRIC OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
4MFA1001	None - 1/8" NPT Female	
4MFA1003	DISS FEMALE w/Hand Tight	
4MFA1006	CHEMETRON (NCG)	
4MFA1005	OHIO/OHMEDA	
4MFA1008	PURITAN BENNETT	

0 - 3 LPM PEDIATRIC OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA3001	None - 1/8" NPT Female	
1MFA3003	DISS FEMALE w/Hand Tight	
1MFA3006	CHEMETRON (NCG)	
1MFA3005	OHIO/OHMEDA	
1MFA3008	PURITAN BENNETT	

0 - 15 LPM OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA1001	None - 1/8" NPT Female	
1MFA1003	DISS FEMALE w/Hand Tight	
1MFA1006	CHEMETRON (NCG)	
1MFA1005	OHIO/OHMEDA	
1MFA1008	PURITAN BENNETT	

0 - 15 LPM OXYGEN FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA1003PTO	DISS FEMALE w/Hand Tight	
1MFA1006PTO	CHEMETRON (NCG)	
1MFA1005PTO	OHIO/OHMEDA	
1MFA1008PTO	PURITAN BENNETT	
1MFA100MPTO	MEDSTAR	

0 - 70 LPM HI FLOW OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA8001	None - 1/8" NPT Female	
1MFA8003	DISS FEMALE w/Hand Tight	
1MFA8006	CHEMETRON (NCG)	
1MFA8005	OHIO/OHMEDA	
1MFA8008	PURITAN BENNETT	

0 - 70 LPM HI FLOW AIR FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER	PRICE
1MFA9001	None - 1/8" NPT Female	
1MFA9003	DISS FEMALE w/Hand Tight	
1MFA9006	CHEMETRON (NCG)	
1MFA9005	OHIO/OHMEDA	
1MFA9008	PURITAN BENNETT	

SELECT FLOWMETER

15 LPM OXYGEN FLOWMETER WITH ADAPTER WITH FLOW SELECTOR BUILT IN

PART NUMBER	ADAPTER	PRICE
3MFA1001	None - 1/8" NPT Female	
3MFA1003	DISS FEMALE w/Hand Tight	
3MFA1006	CHEMETRON (NCG)	
3MFA1005	OHIO/OHMEDA	
3MFA1008	PURITAN BENNETT	

FOR MORE FLOWMETERS, PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM


FLOW SELECTOR



Flow Selector Features

- Give nebulizer treatments without disconnecting the humidifier
- Reduces the changeover time to seconds
- No more loose tubing laying around
- Eliminates tubing nipple adaptors (It has 2 of its own.)
- Fits all flowmeters that have DISS male outlet fittings
- Lightweight

Part # 1000 **Price**



INCREMENT SCALES

0-200cc, 0-1LPM, 0-3 LPM,
0-70 LPM, 0-15 LPM

PRECISION MEDICAL COMPACT 15 LPM FLOWMETERS

OXYGEN



8MFA1005

15 LPM OXYGEN FLOWMETER WITH ADAPTER		
Part #	Adapter	Price
8MFA1001	None - 1/8" NPT Female	
8MFA1002	Diss Female	
8MFA1003	Diss Female w/hand nut	
8MFA1006	Chemetron (NCG)	
8MFA1005	Ohio (Ohmeda)	
8MFA1008	Puritan	
8MFA100M	Medstar	
8MFA1007	Oxequip	



Diss Female



Diss Female with hand nut



8MFA1006PTO

15 LPM OXYGEN FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER		
Part #	Adapter	Price
8MFA1001 PT0	None - 1/8" NPT Female	
8MFA1002 PT0	Diss Female	
8MFA1003 PT0	Diss Female w/hand nut	
8MFA1006 PT0	Chemetron (NCG)	
8MFA1005 PT0	Ohio (Ohmeda)	
8MFA1008 PT0	Puritan	
8MFA100M PT0	Medstar	
8MFA1007 PT0	Oxequip	



Chemetron (NCG)



Y8MFA1006

15 LPM DUAL OXYGEN FLOWMETER WITH ADAPTER		
Part #	Adapter	Price
Y8MFA1001	None - 1/8" NPT Female	
Y8MFA1002	Diss Female	
Y8MFA1003	Diss Female w/hand nut	
Y8MFA1006	Chemetron (NCG)	
Y8MFA1005	Ohio (Ohmeda)	
Y8MFA1008	Puritan	
Y8MFA100M	Medstar	
Y8MFA1007	Oxequip	



Ohio (Ohmeda)



Puritan



Y8MFA1005PT2

15 LPM DUAL OXYGEN FLOWMETER WITH (2) 50 PSI AUXILIARY VALVES AND ADAPTER		
Part #	Adapter	Price
Y8MFA1001 PT2	None - 1/8" NPT Female	
Y8MFA1002 PT2	Diss Female	
Y8MFA1003 PT2	Diss Female w/hand nut	
Y8MFA1006 PT2	Chemetron (NCG)	
Y8MFA1005 PT2	Ohio (Ohmeda)	
Y8MFA1008 PT2	Puritan	
Y8MFA100M PT2	Medstar	
Y8MFA1007 PT2	Oxequip	



Medstar



Oxequip

CHAIN-RING ASSEMBLY



The Chain-Ring Assembly is used to permanently attach Flowmeters to the wall outlet.

Part #	Description	Price
8004	6" Chain w/Closed Ring	
8005	6" Chain w/Open Ring	



50 PSI Auxiliary Valve

PRECISION MEDICAL COMPACT 15 LPM FLOWMETERS

AIR



8MFA2005

15 LPM AIR FLOWMETER WITH ADAPTER		
Part #	Adapter	Price
8MFA2001	None - 1/8" NPT Female	
8MFA2002	Diss Female	
8MFA2003	Diss Female w/hand nut	
8MFA2006	Chemetron (NCG)	
8MFA2005	Ohio (Ohmeda)	
8MFA2008	Puritan	
8MFA200M	Medstar	
8MFA2007	Oxequip	



Diss Female



Diss Female with hand nut



Chemetron (NCG)



Ohio (Ohmeda)



Puritan



Medstar



Oxequip



50 PSI Auxiliary Valve



8MFA2006PTO

15 LPM AIR FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER		
Part #	Adapter	Price
8MFA2001 PT0	None - 1/8" NPT Female	
8MFA2002 PT0	Diss Female	
8MFA2003 PT0	Diss Female w/hand nut	
8MFA2006 PT0	Chemetron (NCG)	
8MFA2005 PT0	Ohio (Ohmeda)	
8MFA2008 PT0	Puritan	
8MFA200M PT0	Medstar	
8MFA2007 PT0	Oxequip	



Y8MFA2005

15 LPM DUAL AIR FLOWMETER WITH ADAPTER		
Part #	Adapter	Price
Y8MFA2001	None - 1/8" NPT Female	
Y8MFA2002	Diss Female	
Y8MFA2003	Diss Female w/hand nut	
Y8MFA2006	Chemetron (NCG)	
Y8MFA2005	Ohio (Ohmeda)	
Y8MFA2008	Puritan	
Y8MFA200M	Medstar	
Y8MFA2007	Oxequip	




Y8MFA2005PT2

15 LPM DUAL AIR FLOWMETER WITH (2) 50 PSI AUXILIARY VALVES AND ADAPTER		
Part #	Adapter	Price
Y8MFA2001 PT2	None - 1/8" NPT Female	
Y8MFA2002 PT2	Diss Female	
Y8MFA2003 PT2	Diss Female w/hand nut	
Y8MFA2006 PT2	Chemetron (NCG)	
Y8MFA2005 PT2	Ohio (Ohmeda)	
Y8MFA2008 PT2	Puritan	
Y8MFA200M PT2	Medstar	
Y8MFA2007 PT2	Oxequip	

FIXED FLOW OXYGEN FLOWMETERS (1/2 - 15 LPM)

Position Increments (LPM)
Off, .50, 1, 2, 3, 4, 5, 6, 8, 10, 15



It's Durable - The lack of a flow tube and the durable metal body prevents breakage problems long associated with standard flowmeters.

It's Compact - The compact size allows the unit to fit into tight headwall areas and makes it virtually damage proof when other equipment is installed.

It's Accurate - The patented multi-spring technology insures accuracy at every flow setting without the need for a visible float as in Thorpe tube flowmeters

It's User-Friendly - The knurled flow adjustment knob and the large, easy-to-read numbers make flow adjustment easy.

It's Versatile - The unit is available for headwall outlets as well as for ceiling mount drops. Both air and oxygen models are available.

FOR MORE FLOWMETERS, PLEASE VISIT
OUR WEBSITE AT:
WWW.RCMEDICAL.COM



	Adapters	Side Connection	Top Connection
	1/8" NPT FE	7MFA1001	7MFA1001TC
	DISS Female Hand Tight	7MFA1003	7MFA1003TC
	Ohmeda	7MFA1005	7MFA1005TC
	Chemtron (NCG)	7MFA1006	7MFA1006TC
	Puritan-Bennett	7MFA1008	7MFA1008TC

NEB BLOCK w/FLOWMETER

Features:

- The Neb Block allows a gas outlet to simultaneously or individually provide various respiration treatments without disconnecting any tubing or equipment.
- Provides easy to use On/Off switch to quickly access oxygen flow to a hand held nebulizer or resuscitation bag.
- Reliable fixed flow orifices provide six ranges of flow for aerosol therapy or resuscitation.



6-7 LPM 8-9 LPM 14-16 LPM
7-8 LPM 10-11 LPM

	6-7 LPM	8MFA6065	
	7-8 LPM	8MFA6075	
	8-9 LPM	8MFA6085	
	10-11 LPM	8MFA6105	
	14-16 LPM	8MFA6155	
	6-7 LPM	8MFA6066	
	7-8 LPM	8MFA6076	
	8-9 LPM	8MFA6086	
	10-11 LPM	8MFA6106	
	14-16 LPM	8MFA6156	
	6-7 LPM	8MFA6063	
	7-8 LPM	8MFA6073	
	8-9 LPM	8MFA6083	
	10-11 LPM	8MFA6103	
	14-16 LPM	8MFA6153	
	6-7 LPM	8MFA6068	
	7-8 LPM	8MFA6078	
	8-9 LPM	8MFA6088	
	10-11 LPM	8MFA6108	
	14-16 LPM	8MFA6158	

SPECIALTY OXYGEN FLOWMETERS



15 LPM O2 Flowmeter for Chemtron Wall Outlet with 50 PSI-Chemtron Female Quick Connect
SFO-5



15 LPM O2 Flowmeter for Chemtron Wall Outlet with 50 PSI-Chemtron Female Quick Connect
SFO-2



Dual 15 LPM Oxygen Flowmeter with Chemtron Quick Connect Inlet and (2) Chemtron Couplers (Auxiliary Valves)
SFO-12



15 LPM O2 Flowmeter for Ohmeda Wall Outlet with 50 PSI-Ohmeda Female Quick Connect
SFO-4 \$



15 LPM O2 Flowmeter for Ohmeda Wall Outlet with 50 PSI-Ohmeda Female Quick Connect
SFO-1



Dual 15 LPM Oxygen Flowmeter with Ohmeda Quick Connect Inlet and (2) Ohmeda Couplers (Auxiliary Valves)
SFO-11

SMALL CYLINDER CARTS (D&E)

#6105.....

41" high

Nonslip Ribbed Hand Grip

Gurney/ Mask Hook

Non-marring Thumb Screw
Front Carrying Handle



#PX-1060.....

OUR BEST SELLING CART

Adjustable height small cylinder cart for "D" and "E" cylinders.



#6103-4.....

41" high
3" casters swivel w/brake



These four-wheeled units offer better control and easy maneuverability in and out of tight spaces. They require no tilting or balancing, giving you total control over the cylinder load.

#1060-F



Holds 4

#6205



Holds 2

LARGE CYLINDER CARTS AND STANDS (H,K)

Dual
6214



Single
6114



All models have 4 wheels.

2 Wheels Only



Part # HCC2-M



610FW-E.....



Fiberglass Cylinder Bases
65045.....

Part #	Description	Price
6114	Single (Most Popular)	
6114-PN	Single (Pneumatic Wheels)	
6214	Dual	
6214-PN	Dual (Pneumatic Wheels)	

SMALL CYLINDER CARTS (D&E)

2-Wheel Cart



Part #	Cylinders Held	Price	Dimensions (Height-Depth-Width)
6041	4		42" - 10" - 18"
6061	6		42" - 15" - 18"
6121	12		42" - 20" - 22"

With Sound Deadening Kit (below)

6041-SDK	4		42" - 10" - 18"
6061-SDK	6		42" - 15" - 18"
6121-SDK	12		42" - 20" - 22"

4-Wheel Cart



Part #	Cylinders Held	Price	Dimensions (Height-Depth-Width)
6064	6		42" - 15" - 18"
6124	12		42" - 20" - 22"
6184*	18		42" - 30" - 16"
6244*	24		42" - 30" - 20"

With Sound Deadening Kit (below)

6124-SDK	12		42" - 20" - 22"
6184-SDK*	18		42" - 30" - 16"
6244-SDK*	24		42" - 30" - 20"

* -4" wheels

Multiple D & E Size Cylinder Stands With Wheels



Part #	Cylinders Held	Price	Dimensions (Height-Depth-Width)
6040-PC	4		23" - 13" - 10"
6060-PC	6		23" - 12" - 15"
6120-PC	12		23" - 15" - 22"

With Sound Deadening Kit (below)

6040-PC-SDK	4		23" - 13" - 11"
6060-PC-SDK	6		23" - 12" - 15"
6120-PC-SDK	12		23" - 15" - 22"

ALL CARTS ABOVE ARE PAINTED GREEN

RED CARTS - ADD A "R" TO END OF PART NUMBER, PRICE IS \$ MORE
YELLOW CARTS - ADD A "Y" TO END OF PART NUMBER, PRICE IS \$ MORE

SOUND DEADENING KITS (for Carts Above)



Plastic Bumper Guards

Sound Deadening Matting

Anthony has developed an exclusive sound deadening system for their multiple cylinder carts, used to distribute cylinders in and around hospitals, rest homes, convalescent homes, etc. or where QUIET is essential.

These sound deadening kits include heavy-duty sound deadening matting and plastic bumper guards. This formula eliminates the rattling & clanging of steel & aluminum cylinders while in transport.

**SOUND DEADENING KITS ALSO SOLD SEPARATELY
 PLEASE VISIT OUR WEBSITE: WWW.RCMEDICAL.COM**



10" x 14" Plastic Sign
 Part # SP-51



10" x 14" Plastic Sign
 Part # SP-52

SMALL CYLINDER STANDS AND HOLDERS (D&E)

Multiple D & E Size Cylinder Stands



Part #	Cylinders Held	Price	Dimensions (Height-Depth-Width)
6040	4		20" - 13" - 11"
6060	6		20" - 12" - 15"
6120	12		20" - 15" - 22"
6180	18		20" - 15" - 32"
6240	24		20" - 20" - 32"
6400	40		20" - 40" - 25"

With Sound Deadening Kit (below)

6040-SDK	4		20" - 13" - 11"
6060-SDK	6		20" - 12" - 15"
6120-SDK	12		20" - 15" - 22"
6180-SDK	18		20" - 15" - 32"
6240-SDK	24		20" - 20" - 32"
6400-SDK	40		20" - 40" - 25"

Can be bolted to floor or free standing.

ALL STANDS ABOVE ARE PAINTED GREEN

RED STANDS - ADD A "R" TO END OF PART NUMBER, PRICE IS \$ MORE

YELLOW STANDS - ADD A "Y" TO END OF PART NUMBER, PRICE IS \$ MORE



MULTIPLE CYLINDER STANDS

CHROME



PX-1025

PX-1025 "D/E" Rack / 6 cylinders 9" x 15" x 10"

PX-1027 "D/E" Rack / 12 cylinders 9" x 15" x 20"

Bolts to floor or free standing.



For any solid surfaces,
crash carts, more...
Part # 610WM



For wheelchairs
Part # 610WC

FOR MORE INFORMATION ON SMALL CYLINDER BASES,
CYLINDER HOLDERS FOR BEDS, IV POLES, M6 CYLINDER CARTS AND STANDS,
AND MORE, PLEASE VISIT OUR WEBSITE AT:

WWW.RCMEDICAL.COM

OXYGEN REGULATORS - LITERS PER MINUTE

Responsive Respiratory Fixed Flow Oxygen Regulator CGA 870 for Small Cylinder (D & E)




Part #120-3040
for small cylinder (CGA 870)

ALL BRASS

0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 settings



Part #21017.....
for small cylinder (870)



Part #120-1061C.....
for large cylinder (CGA 540)

0, 1, 2, 4, 6, 8, 10, 12, 15, 18, 20, 25 settings



Part #21025.....
1 - 25 LPM

Two High Pressure Outlets Also

Western Medica Oxygen Regulators

(Diaphragm Type)

6 Year Warranty



15 LPM

Single stage regulator for large cylinders.
15 liters per minute.
ORDER #M1-540-15 FG

With Gauge-Savers
ORDER #M1-540-15 FGP



M1-540-15 FM

Single stage flowmeter/regulator for large cylinders.
15 liters per minute.

ORDER #M1-540-15 FM.....

Single stage flowmeter/regulator for small cylinders.
15 liters per minute.

ORDER #M1-870-15 FM.....

 **MADE IN USA** 



15 LPM

Single stage regulator for small cylinders.
15 liters per minute.
ORDER #M1-870-15 FG

With Gauge-Savers
ORDER #M1-870-15 FGP



15 LPM

Single stage, swivel yoke regulator for small cylinders.
15 liters per minute.
ORDER #MR-870-15 FG

With Gauge-Savers
ORDER #MR-870-15 FGP

To order 8 or 5 liter regulators, add an 8 or a 5 to the end of our part number. The same price as our 15 liter regulators.

HIGH PRESSURE REGULATORS



Adjustable 0 - 100 psi with DISS Outlet

Gas	Part #	Cylinder Size	Price
Oxy.	MI 540 PG	CGA 540 (Large)	
Oxy.	MI 870 PG	CGA 870 (Small)	
Air	MI 346 PG	CGA 346 (Large)	
Air	MI 950 PG	CGA 950 (Small)	



Preset 50 psi with DISS Outlet


Gas	Part #	Cylinder Size	Price
Oxy.	MIA 540 P	CGA 540 (Large)	
Air	MIA 346 P	CGA 346 (Large)	
N ₂ O	MIA 326 P	CGA 326 (Large)	
CO ₂	MIA 320 P	CGA 320 (Large)	



Preset 50 psi with DISS Outlet

Gas	Part #	Cylinder Size	Price
Oxy.	MIA 870 P	CGA 870 (Small D&E)	
Air	MIA 950 P	CGA 950 (Small D&E)	
N ₂ O	MIA 910 P	CGA 910 (Small D&E)	
CO ₂	MIA 940 P	CGA 940 (Small D&E)	

MADE IN USA



Preset 50 psi with Ohmeda Quick Connect

Gas	Part #	Cylinder Size	Price
Oxy.	SRO-19	CGA 540 (Large)	
Air	SRO-34	CGA 346 (Large)	



Preset 50 psi with Chemetron Quick Connect

Gas	Part #	Cylinder Size	Price
Oxy.	SRO-21	CGA 540 (Large)	
Air	SRO-36	CGA 346 (Large)	



Preset 50 psi with Ohmeda Quick Connect


Gas	Part #	Cylinder Size	Price
Oxy.	SRO-30	CGA 870 (Small D&E)	
Air	SRO-32	CGA 950 (Small D&E)	



Preset 50 psi with Chemetron Quick Connect

Gas	Part #	Cylinder Size	Price
Oxy.	SRO-29	CGA 870 (Small D&E)	
Air	SRO-31	CGA 950 (Small D&E)	


SPECIALTY OXYGEN REGULATORS



Fixed flow Regulator
With Ohmeda 50 psi
Outlet. .5 -25 lpm
0 - 25 LPM
Part # SRO-11.....



Fixed flow Regulator
With Chemetron 50 psi
Outlet. .5 - 25 lpm
Part # SRO-12.....



Fixed flow Regulator
With Puritan-Bennett
50 psi outlet. .5 - 25 lpm
Part # SRO-13.....



Oxygen Regulator/
Flowmeter 15 LPM,
Large Cylinder with
50 PSI Valve
Part # SRO-14.....



Oxygen Regulator with
Dual Chemetron Female
Couplers, 50 PSI, Large
Cylinder
Part # SRO-15.....



Oxygen Regulator with
Dual Ohmeda Female
Couplers, 50 PSI, Large
Cylinder
Part # SRO-16.....



Large Cylinder Oxygen
Regulator with (2)
15 LPM Flowmeters
Part # SRO-17.....



Large Cylinder Oxygen
Regulator with (2)
15 LPM Flowmeters and
(2) Auxiliary Valves
Part # SRO-18.....



Oxygen Regulator, Preset
50 PSI with Ohmeda
Quick Connect, Large
Cylinder
Part # SRO-19.....



Oxygen Regulator, Preset
50 PSI with Puritan
Bennett Quick Connect,
Large Cylinder
Part # SRO-20.....



Oxygen Regulator, Preset
50 PSI with Chemetron
Quick Connect, Large
Cylinder
Part # SRO-21.....



Oxygen Regulator with
15 LPM flowmeter and
50 PSI Ohmeda Quick
Connect, Large Cylinder
Part # SRO-22.....



Oxygen Regulator with
15 LPM Flowmeter and
50 PSI Chemetron Quick
Connect, Large Cylinder
Part # SRO-23.....



Oxygen Regulator, Preset
50 PSI with Chemetron
Quick Connect, Small
Cylinder (CGA 870)
Part # SRO-29.....



0 - 25 LPM Regulator
with Chemetron Female
50 PSI Outlet for Large
Cylinder (CGA 540)
Part # SRO-25.....




0 - 25 LPM Regulator
with Ohmeda Female
50 PSI Outlet for Large
Cylinder (CGA 540)
Part # SRO-26.....

Specialty AIR Regulators, Visit our Website at WWW.RCMEDICAL.COM

For Regulator Accessories See Page 31

HELIOX, AIR, CO₂ REGULATORS

Heliox Regulators



0 - 20 LPM


Part #	Ratio	For Cylinder size	Price
2210-20	70 - 30	Small D and E CGA 890 Yoke	
2211-20	80 - 20	Small D and E CGA 890 Yoke	
2212-20	70 - 30	Large CGA 280 Nut/Nipple	
2213-20	80 - 20	Large CGA 280 Nut/Nipple	

WESTERN MEDICA MEDICAL REGULATORS (AIR or CARBON DIOXIDE)


Single-stage regulators preset to 50 psi or adjustable to 100 psi, or with a flowmeter. A sensitive neoprene diaphragm assembly provides optimum accuracy and sensitivity. Dual inlet filtering. Available for small or large cylinders.

		AIR	CO ₂
Preset 50 PSI		M1A-346P for large cylinder	M1-320P for large cylinder
Adjustable 0 - 100 PSI		M1-346PG for large cylinder	M1-320PG for large cylinder
12 Liters Per Minute (lpm)		M1-346-15FM for large cylinder	M1-320-12FM for large cylinder
Preset 50 PSI		M1A-950P for small cylinder	M1-940P for small cylinder
Adjustable 0 - 100 PSI		M1-950PG for small cylinder	M1-940PG for small cylinder
12 Liters Per Minute (lpm)		M1-950-15FM for small cylinder	M1-940-12FM for small cylinder

EASTERN RAIL SYSTEMS



Simple Outlet Expansion allows your hospital to use its existing outlets, and to create the number and placement of outlets needed. Simply mount a length of our rail on the wall - attach the needed outlets to the rail with our thumb screw clamp - attach a length of high-pressure hose to the block, and run it back to your existing outlet.




- Easy external mounting - no need to shut down the unit.
- No more weight on outlets
- Shift the weight of humidifiers & suction bottles to the rail
- Eliminate overloading
- Reduce outlet leakage. Save time and money
- Increase efficiency
- Works with all piping systems
- Have outlets where you need them NOW!


Kit comes with:

- 1 each 2 Foot Rail
- 2 each Mounting Hardware
- 1 each 3 Foot Hose with Wall Adapter
- 1 each Double Gas Block

Outlet connection	GAS	Part #	Price
DISS	Oxygen	ERO-1DS	
Chemetron	Oxygen	ERO-1CH	
Ohmeda	Oxygen	ERO-10M	
DISS	AIR	ERA-1DS	
Chemetron	AIR	ERA-1CH	
Ohmeda	AIR	ERA-10M	
DISS	VAC	ERV-1DS	
Chemetron	VAC	ERV-1CH	
Ohmeda	VAC	ERV-10M	



**FOR MORE EASTERN RAIL KITS
AND ACCESSORIES,
PLEASE GO TO OUR WEBSITE AT:
WWW.RCMEDICAL.COM**



CUSTOM HOSE ASSEMBLIES

Ordering is easy as 1, 2, 3 . . . 4


How to Order:

Western Medica manufactures custom hose assemblies that conform to CGS, and NFPA standards. All assemblies are pressure-tested. Select your requirement from the charts below.

Example: OXY-024-FNN-OHF

- 1) OXY Oxygen Service
- 2) 024 24" Long
- 3) FNN Inlet Connection
- 4) OHF Outlet Connection

Evacuation Hose (purple) price may be more, depending on your adaptor



(Pricing for the above example)

FNN
OHF
Hose (.35 per inch x 24")
Assembly Charge

Gas Service

Gas	Service	Color
OXY	Oxygen	Green
AIR	Medical Air	Yellow
SUC	Suction	White
N2O	Nitrous Oxide	Blue
EVC	Evacuation	Purple
CO ₂	Carbon Dioxide	Gray
NIT	Nitrogen	Black
GNR	General	Clear

Adapters pictured below may be ordered Separately.

When ordering adapters only, please specify gas.
Example: chm - oxy

MINIMUM HOSE LENGTH — 1 FOOT

HOSE ORDERS FOR CUSTOM HOSES MAY NOT BE CANCELLED

All hose assemblies are 1/4" conductive (unless specified)

#CHM Chemetron (NCG) Male	#FNN Diss Female with Hex nut	#FHT Diss Female with Hand Nut	#OHM Ohmeda Male
#MBA Diss Male	#SHM Schrader Male	#OXM Oxequip Male	#PBM Puritan - Bennett Male
#CHF Chemetron (NCG) Female	#SHF Schrader Female	#PBF Puritan - Bennett Female	#OHF Ohmeda Female
#DBA Demand Body Adapter	#MP2 Male Pipe 1/8"	#MP4 Male Pipe 1/4"	
#FP2 Female Pipe 1/8"	#FP4 Female Pipe 1/4"	#BLK Blank (No Connection)	

REINFORCED MEDICAL HOSE



Hose only! per foot

Conductive, Non-UL Listed
400 psi Working Pressure
1500 psi Burst Pressure

MH4-04C	1/4" ID, Nitrous Oxide (Blue)
MH4-08C	1/4" ID, Carbon Dioxide (Grey)
MH4-12C	1/4" ID, Nitrogen (Black)
MH4-16C	1/4" ID, Medical Air (Yellow)
MH4-22C	1/4" ID, Suction (White)
MH4-24C	1/4" ID, Oxygen (Green)
MH5-22C	5/16" ID, Suction (White)

FERRULE



For 1/4" hose
.525 ID
#M-625-A.....

VICE-GRIPS

For 3/16", 1/4" and
5/16" ID Hose
C-5A



COMMON HOSE ASSEMBLIES

DISS FEMALE HANDTIGHT X DISS FEMALE HANDTIGHT



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-100	
5 Foot Air	CH-101	
10 Foot Oxygen	CH-102	
10 Foot Air	CH-103	
15 Foot Oxygen	CH-104	
15 Foot Air	CH-105	
20 Foot Oxygen	CH-106	
20 Foot Air	CH-107	

DISS FEMALE HANDTIGHT X OHMEDA MALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-132	
5 Foot Air	CH-133	
10 Foot Oxygen	CH-134	
10 Foot Air	CH-135	
15 Foot Oxygen	CH-136	
15 Foot Air	CH-137	
20 Foot Oxygen	CH-138	
20 Foot Air	CH-139	

DISS FEMALE HANDTIGHT X CHEMETRON MALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-108	
5 Foot Air	CH-109	
10 Foot Oxygen	CH-110	
10 Foot Air	CH-111	
15 Foot Oxygen	CH-112	
15 Foot Air	CH-113	
20 Foot Oxygen	CH-114	
20 Foot Air	CH-115	

DISS FEMALE HANDTIGHT X OHMEDA FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-140	
5 Foot Air	CH-141	
10 Foot Oxygen	CH-142	
10 Foot Air	CH-143	
15 Foot Oxygen	CH-144	
15 Foot Air	CH-145	
20 Foot Oxygen	CH-146	
20 Foot Air	CH-147	

DISS FEMALE HANDTIGHT X CHEMETRON FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-116	
5 Foot Air	CH-117	
10 Foot Oxygen	CH-118	
10 Foot Air	CH-119	
15 Foot Oxygen	CH-120	
15 Foot Air	CH-121	
20 Foot Oxygen	CH-122	
20 Foot Air	CH-123	

OHMEDA MALE X OHMEDA FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-148	
5 Foot Air	CH-149	
10 Foot Oxygen	CH-150	
10 Foot Air	CH-151	
15 Foot Oxygen	CH-152	
15 Foot Air	CH-153	
20 Foot Oxygen	CH-154	
20 Foot Air	CH-155	

CHEMETRON MALE X CHEMETRON FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-124	
5 Foot Air	CH-125	
10 Foot Oxygen	CH-126	
10 Foot Air	CH-127	
15 Foot Oxygen	CH-128	
15 Foot Air	CH-129	
20 Foot Oxygen	CH-130	
20 Foot Air	CH-131	

PURITAN BENNETT MALE X PURITAN BENNETT FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-156	
5 Foot Air	CH-157	
10 Foot Oxygen	CH-158	
10 Foot Air	CH-159	
15 Foot Oxygen	CH-160	
15 Foot Air	CH-161	
20 Foot Oxygen	CH-162	
20 Foot Air	CH-163	

DISS FEMALE HANDTIGHT X PURITAN BENNETT MALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-164	
5 Foot Air	CH-165	
10 Foot Oxygen	CH-166	
10 Foot Air	CH-167	
15 Foot Oxygen	CH-168	
15 Foot Air	CH-169	
20 Foot Oxygen	CH-170	
20 Foot Air	CH-171	

DISS FEMALE HANDTIGHT X PURITAN BENNETT FEMALE



DESCRIPTION	PART #	PRICE
5-Foot Oxygen	CH-172	
5 Foot Air	CH-173	
10 Foot Oxygen	CH-174	
10 Foot Air	CH-175	
15 Foot Oxygen	CH-176	
15 Foot Air	CH-177	
20 Foot Oxygen	CH-178	
20 Foot Air	CH-179	

HUDSON RCI DISPOSABLE ADAPTERS AND CONNECTORS



No. 1667.....
Ventilator Monitoring Adaptor Circuit

- Allows monitoring of ventilatory weaning parameters in intubated patients
- Includes oxygen nipple adaptor, aerosol tee, 6" CORR-A-FLEX tube, 15mm/22mm adaptor, m/m one-way valve and m/f one-way valve
- 50 per case



No. 1637.....
Thermometer with Tee Adapter

- Combines No. 1635 Thermometer and No. 1077 Aerosol Tee Connector
- Packaged 50 to a case



No. 1743.....
Valved Neb Tee

- 22mm O.D. x 18mm I.D. x 22mm I.D. connections
- Tee may be left in-line when nebulizer is removed
- Port accepts most small volume nebulizer
- 30 per case



No. 1648.....
Aerosol Tee Connector

- 22mm O.D. x 15mm I.D./22mm O.D. x 22mm O.D.
- 50 per case



No. 1641.....
Elbow Adaptor

- Standard 22mm I.D. connector on vertical axis
- 22mm O.D. connection on horizontal axis
- Packaged 50 to a case



No. 1642.....
Pressure Line Adaptor

- 22mm I.D./O.D. connections
- Barbed elbow accepts 5 to 7mm I.D. tubing and swivels 360 degrees
- Packaged 50 to a case



No. 1623.....
Patient Elbow Adaptor

- 15mm O.D. x 15mm I.D./22mm O.D.
- 50 per case



No. 1644.....
One-Way Valve

- Standard 22 mm I.D./O.D. connections
- Valve offers low flow resistance and a capped monitoring port
- Packaged 50 to a case




No. 1638.....
Aerosol Tee Connector

- Horizontal ports of the tee are 22mm O.D.
- Vertical port is 18mm I.D./22mm O.D.
- Packaged 50 to a case



No. 1077.....
Aerosol Tee Connector

- Anti-spill design minimizes aspiration of liquids
- Horizontal axis has standard 22mm O.D. connections
- Vertical axis incorporates a standard 15mm I.D. connection and two 22mm O.D. concentric rings
- Packaged 50 to a case



No. 1639.....
Nebulizer Tee Connector

- Standard 22 mm I.D./O.D. connections
- Nebulizer port is 18mm I.D. which accepts most small volume nebulizers
- Packaged 50 to a case



No. 1421.....
Universal Cuff Adaptor

- Standard 22mm I.D. connections
- Packaged 50 to a case



No. 1422.....
Multi Adaptor

- Standard 15mm I.D./22mm O.D. connections
- Packaged 50 to a case



No. 1409.....
Large Tubing Connector

- Standard 22mm O.D. connections
- Packaged 50 to a case



No. 1640.....
Tubing Adaptor

- Standard 22mm I.D./O.D.
- Packaged 50 to a case



No. 1751.....
MDI Adaptor

- 15mm O.D x 15mm I.D. connections for in-line placement in 15mm tubing
- 50 per case



No. 1659.....
MDI Adaptor

- 22mm O.D x 22mm I.D. connections for in-line placement in 22mm tubing
- 50 per case



No. 1423.....
Tubing Adaptor

- 22mm I.D. x 5 - 7mm tapered nipple
- Packaged 50 to a case



No. 1669.....
Monitoring Adaptor

- 15mm O.D. x 15mm I.D./22mm O.D.
- Standard luer fitting
- 50 per case

DISPOSABLES

#2200 Flexible Connector (Silicone Rubber)



22mm I.D. x 22mm I.D.

This adapter was designed almost identical to the ones listed below:

INTEC 2200 MARQUEST 172200
 BENNETT 3443 AM SYSTEMS 1139
 Instrumentation Ind. 101-00-3443

Part Number 2200
per case of 50
each

27mm I.D.



22mm I.D.

Tapered Silicone Rubber Connector

27mm ID x 22mm ID

Part Number 1426
per case of 10



MDI Accessory by Hudson

The POCKET CHAMBER is equipped with a highly sensitive silicone valve that opens easily without restricting the air flow. This one-way valve helps keep patients from accidentally blowing medication out of the chamber, and may aid patients experiencing difficulty timing their inhalation.

Part Number 1001-10
per case of 10



Main Flow Bacterial / Viral Filter

- 22mm ID x 15mm ID/22mm OD
- Bacterial filtration efficiency 99.999+%
- Viral filtration efficiency 99.99+%
 - Dead space: 42mL
- Flow resistance: approximately 1.5 cmH₂O at 60 LPM
- Individually packaged 50 to a case

Part Number 1605
per case of 50



Supply Line Adapter

For extending smooth bore oxygen tubing
 Tapered for excellent gripping

Part Number 65100
per case of 200



Corrugated Tubing

- Continuous 100-foot roll of tubing with cuttable 6-inch sections
- Packaged in protective, convenient, plastic-lined dispenser box

Part Number 1680
per roll



Ventilator Circuits

Adult Standard Ventilator Circuit, 72" long.
 Use on Puritan-Bennett 7200A; Servo 300, 900A, 900B, and 900C; Hamilton Veolar; Emerson; Engstrom Erica; PPG Irisa; Infrasonics Adult Star and Bird 8400ST

Part Number 1613
per case of 10



Voldyne Incentive Spirometer

Voldyne volumetric exerciser features Good - Better - Best flow window in a compact, ergonomic design with advanced, low work-of-breathing filter design and built in handle. 5000 ml

Part Number 8884719009
per case of 12



E-Z Wrap

Soft closed-cell foam tubes designed to prevent ear chafing.
 Light gray color
 Packed 50 pair per case.

Part Number 1016
per case of 50



Noseclip

Part Number 64019
per case of 100



Aerosol Drainage Bags

Disposable Aerosol Drainage Bags with tee adaptor fits 22mm ID corrugated tubing.
 Bag has plastic chain for securing and convenient drainage port. 750 cc capacity.
 With wye adaptor.

Part Number 1742
per case of 50

RESPIROMETERS

MARK 8



1.4" Face

MARK 14



1.4" Face

HALOSCALE STANDARD



2.25" Face

For measurement of minute volumes as well as tidal volumes. Two dials. Small dial records 0-1 liter. Large dial 0-100 liters.

Specially suited for the measurement of tidal volumes and may also be used for measurement of minute volumes. Two dials. Small dial records 0-10 liters. Large dial records 0-1 liter.

The original design of Haloscale Respirometer incorporating concentric dials and an enlarged face for improved readability. Equally convenient for measuring both tidal and minute volumes. Large 57 mm dial diameter. Red/black hands and black scales on a white background with yellow highlight ring. Two concentric outer dials. First dial records 0-1 liter. Second dial records 0-100 liters.

ORDER #	RESPIROMETERS DESCRIPTION	PRICE
2750.....	Mark 8 Respirometer	
2751.....	Mark 14 Respirometer	
2768.....	Haloscale Standard.....	

RESPIROMETERS WITH SPIRO GUARD

2756.....	Mark 8 Respirometer	
2757.....	Mark 14 Respirometer	
2769.....	Haloscale Standard.....	

RESPIROMETERS WITH SPIRO GUARD AND EXTRA LARGE CARRYING CASE

2762.....	Mark 8 Respirometer	
2763.....	Mark 14 Respirometer	
2772.....	Haloscale Standard.....	

Used Respirometers

Fully Serviced, Cleaned,
Oiled and Calibrated

1 Year Warranty

Order #2799
(Case included)

ELECTRONIC DIGITAL RESPIROMETER

A respirometer long overdue in this modern electronic age



ANESTHESIA ASSOCIATES

MODEL 295 DIGITAL RESPIROMETER IS AN ELECTRONIC VERSION OF THE MECHANICAL RESPIROMETER. IT DOES THE SAME JOB PLUS

- Easier to read
- Saving on the initial purchase and maintenance
- More reliable
- Extremely durable
- No gears, shafts, or bearings to break down in the electronic unit
- Separates for ease of patient use, sterilization and cleaning

Part #		Price
295	Respirometer	
Options		
296	6' Retractable Cable (extended length)	
298	Carrying Case	
299	Spare Flow Sensing Unit	

Model 295, a two piece detachable respirometer comes complete ready to use with 9v battery and one (1) 22 mm female x 22 mm female adapter.



Rubber Protective Cover for Small Faced Respirometers (mk. 8, 14)

















Part # 2700 Spiro-Guard \$

REPAIRS





Includes cleaning, calibration and oiling. Parts are extra

All Models Repaired including:
Wright • Ferraris • Haloscale
BOC • Fraser - Harlake • Ohmeda
Boehringer







NPT PIPE FITTINGS

<p>Female Tee</p>  <p>Part No. 1/8" NPT..... 7034 1/4" NPT..... 7035</p>	<p>Street Tee</p>  <p>Part No. 1/8" NPT..... 7032 1/4" NPT..... 7033</p>	<p>Cross</p>  <p>Part No. 1/8" NPT..... 7036 1/4" NPT..... 7037</p>	<p>45° Street Elbow</p>  <p>Part No. 1/8" NPT..... 7056 1/4" NPT..... 7057</p>
<p>Connectors</p>  <p>Part No. 1/8" Female x 1/8" Male 1017 1/8" Female x 1/4" Male 7024 1/4" Female x 1/8" Male..... 7025 1/4" Female x 1/4" Male 7026</p>	<p>Hex Nipple</p>  <p>Part No. 1/8" Male x 1/8" Male..... 7020 1/4" Male x 1/4" Male..... 7021 1/8" Male x 1/4" Male..... 7022</p>	<p>Hose Barb Connector</p>  <p>Part No. 1/8" Male x 1/8" Dia. 7038 1/8" Male x 3/16" Dia..... 7039 1/8" Male x 1/4" Dia. 7040 1/4" Male x 1/4" Dia. 7043</p>	<p>Coupler</p>  <p>Part No. 1/8" Female x 1/8" Female 7027 1/4" Female x 1/4" Female 7028 1/4" Female x 1/8" Female 7029</p>
<p>90° Elbow</p>  <p>Part No. 1/8" Female x 1/8" Female 7058 1/4" Female x 1/4" Female 7059</p>	<p>Plug</p>  <p>Part No. 1/4" NPT..... 7017 1/8" NPT..... 7019</p>	<p>Hose Barb Connector</p>  <p>Part No. 1/8" Female x 1/4" Dia. 7046 1/4" Female x 1/4" Dia. 7049</p>	<p>Cap</p>  <p>Part No. 1/4" NPT..... 7016 1/8" NPT..... 7018</p>
<p>Pipe Tee</p>  <p>Part No. 1/8" NPT..... 7030 1/4" NPT..... 7031</p>	<p>Drop Ear</p>  <p>Part No. Tee 1/8" F x 1/8" F x 1/8" F 7081 Elbow 1/8" F x 1/8" F..... 7082</p>	<p>Pipe Nipple</p>  <p>Part No. 1/8" Male x 1 1/2" Long 7050 1/8" Male x 2" Long 7051 1/8" Male x 2 1/2" Long 7052 1/4" Male x 1 1/2" Long 7053 1/4" Male x 2" Long 7054 1/4" Male x 2 1/2" Long 7055 1/8" Male x 3 1/2" Long 7062</p>	<p>90° Elbow</p>  <p>Part No. 1/8" F x 1/8" M..... 7063 1/4" F x 1/4" M..... 7065</p>



"Y" BLOCK ASSEMBLIES

Ohmeda Style (Ohio)	Chemetron Style (NCG)	Puritan-Bennett Style	Ohmeda Style "Y" Assembly--"Stubby"																																				
																																							
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td>2152</td><td>O₂</td><td></td></tr> <tr><td>2252</td><td>Air</td><td></td></tr> <tr><td>2352</td><td>Vac</td><td></td></tr> </table>	2152	O ₂		2252	Air		2352	Vac		<table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td>3152</td><td>O₂</td><td></td></tr> <tr><td>3252</td><td>Air</td><td></td></tr> <tr><td>3352</td><td>Vac</td><td></td></tr> </table>	3152	O ₂		3252	Air		3352	Vac		<table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td>4152</td><td>O₂</td><td></td></tr> <tr><td>4252</td><td>Air</td><td></td></tr> <tr><td>4352</td><td>Vac</td><td></td></tr> </table>	4152	O ₂		4252	Air		4352	Vac		<table border="1" style="width: 100%; border-collapse: collapse;"> <tr><td>2162</td><td>O₂</td><td></td></tr> <tr><td>2262</td><td>Air</td><td></td></tr> <tr><td>2362</td><td>Vac</td><td></td></tr> </table>	2162	O ₂		2262	Air		2362	Vac	
2152	O ₂																																						
2252	Air																																						
2352	Vac																																						
3152	O ₂																																						
3252	Air																																						
3352	Vac																																						
4152	O ₂																																						
4252	Air																																						
4352	Vac																																						
2162	O ₂																																						
2262	Air																																						
2362	Vac																																						


T BLOCKS

OHMEDA Style											
	Gas	Part #	Price		Gas	Part #	Price		Gas	Part #	Price
	Oxy	2152-T			Oxy	2152-TL			Oxy	2152-TR	
	Air	2252-T			Air	2252-TL			Air	2252-TR	
	Vac	2352-T			Vac	2352-TL			Vac	2352-TR	
CHEMETRON Style											
	Oxy	3152-T			Oxy	3152-TL			Oxy	3152-TR	
	Air	3252-T			Air	3252-TL			Air	3252-TR	
	Vac	3352-T			Vac	3352-TL			Vac	3352-TR	

TRI-BLOCK ASSEMBLIES

	Ohmeda Style		
	SP-112	O ₂	
	SP-113	Air	
	SP-114	Vac	
	Chemetron Style		
	SP-115	O ₂	
	SP-116	Air	
	SP-117	Vac	

DISS "Y" ASSEMBLY WITH CHECK VALVES

	7211	O ₂	
(7211)	7212	Air	
	7213	Vac	
<p style="margin: 0;">FOR MORE "Y" BLOCK ASSEMBLIES PLEASE VISIT OUR WEBSITE WWW.RCMEDICAL.COM</p>			

DISS "Y" WITH 2 OHMEDA COUPLERS

	SP-214	O ₂	
	SP-214A	Air	
	SP-214V	Vac	

DISS "Y" WITH 2 CHEMETRON COUPLERS

	SP-215	O ₂	
	SP-215A	Air	
	SP-215V	Vac	

QUICK CONNECTS (MALE)

Ohmeda

1/8" NPT Male



Gas	Part #	Price
Oxygen.....	2101.....	
Air.....	2201.....	
Vacuum.....	2301.....	

1/8" NPT Female



Oxygen.....	2102.....	
Air.....	2202.....	
Vacuum.....	2302.....	

1/4" NPT Male



Oxygen.....	2103.....	
Air.....	2203.....	
Vacuum.....	2303.....	

DISS Male



Oxygen.....	2105.....	
Air.....	2205.....	
Vacuum.....	2305.....	

DISS Hex Nut



Oxygen.....	2106.....	
Air.....	2206.....	
Vacuum.....	2306.....	

DISS Hand Tight



Oxygen.....	2108.....	
Air.....	2208.....	
Vacuum.....	2308.....	

1/4" Hose Barb



Oxygen.....	2111.....	
Air.....	2211.....	
Vacuum.....	2311.....	

DISS Male with Check Valve



Oxygen.....	2114.....	
Air.....	2214.....	
Vacuum.....	2314.....	4

Chemetron

1/8" NPT Male



Gas	Part #	Price
Oxygen.....	3101.....	
Air.....	3201.....	
Vacuum.....	3301.....	

1/8" NPT Female



Oxygen.....	3102.....	
Air.....	3202.....	
Vacuum.....	3302.....	

1/4" NPT Male



Oxygen.....	3103.....	
Air.....	3203.....	
Vacuum.....	3303.....	

DISS Male



Oxygen.....	3105.....	
Air.....	3205.....	
Vacuum.....	3305.....	

DISS Hex Nut



Oxygen.....	3106.....	
Air.....	3206.....	
Vacuum.....	3306.....	

DISS Hand Tight



Oxygen.....	3108.....	
Air.....	3208.....	
Vacuum.....	3308.....	

1/4" Hose Barb



Oxygen.....	3111.....	
Air.....	3211.....	
Vacuum.....	3311.....	

DISS Male with Check Valve



Oxygen.....	3114.....	
Air.....	3214.....	
Vacuum.....	3314.....	

Puritan Bennett

1/8" NPT Male



Gas	Part #	Price
Oxygen.....	4101.....	
Air.....	4201.....	
Vacuum.....	4301.....	

1/8" NPT Female



Oxygen.....	4102.....	
Air.....	4202.....	
Vacuum.....	4302.....	

1/4" NPT Male



Oxygen.....	4103.....	
Air.....	4203.....	
Vacuum.....	4303.....	

DISS Male with Knob



Oxygen.....	4105.....	
Air.....	4205.....	
Vacuum.....	4305.....	

DISS Hex Nut



Gas	Part #	Price
Oxygen.....	4106.....	
Air.....	4206.....	
Vacuum.....	4306.....	

DISS Hand Tight



Oxygen.....	4108.....	
Air.....	4208.....	
Vacuum.....	4308.....	

1/4" Hose Barb w/Knob



Oxygen.....	4111.....	
Air.....	4211.....	
Vacuum.....	4311.....	

DISS Male w/Knob & Check Valve



Oxygen.....	4114.....	
Air.....	4214.....	
Vacuum.....	4314.....	

QUICK CONNECTS (FEMALE)

CHEMETRON STYLE CHECK UNITS (Female Couplers)

1/8" NPT Female



Gas	Part No.
O ₂	3150
Air.....	3250
Vac*	3350

DISS Male



Gas	Part No.
O ₂	3155
Air.....	3255
Vac*	3355

1/4" NPT Male



Gas	Part No.
O ₂	3153
Air.....	3253
Vac*	3353

DISS Female



Gas	Part No.
Hand Tight O ₂	3158
Hand Tight Air.....	3258
Hand Tight Vac*	3358


OHMEDA STYLE CHECK UNITS (Female Couplers)

1/8" NPT Female



Gas	Part No.
O ₂	2150
Air.....	2250
Vac*	2350

DISS Male



Gas	Part No.
O ₂	2155
Air.....	2255
Vac*	2355

1/4" NPT Male



Gas	Part No.
O ₂	2153
Air.....	2253
Vac*	2353

DISS Female



Gas	Part No.
Hand Tight O ₂	2158
Hand Tight Air.....	2258
Hand Tight Vac*	2358

PURITAN-BENNETT STYLE CHECK UNITS (Female Couplers)

1/8" NPT Female



Gas	Part No.
O ₂	4150
Air.....	4250
Vac*	4350

DISS Male



Gas	Part No.
O ₂	4155
Air.....	4255
Vac*	4355

1/4" NPT Male



Gas	Part No.
O ₂	4153
Air.....	4253
Vac*	4353

DISS Female



Gas	Part No.
Hand Tight O ₂	4158
Hand Tight Air.....	4258
Hand Tight Vac*	4358

SCHRADER STYLE CHECK UNITS (Female Couplers)

1/8" NPT Female



Gas	Part No.
O ₂	6150
Air.....	6250
Vac*	6350

DISS Male



Gas	Part No.
O ₂	6155
Air.....	6255
Vac*	6355

1/4" NPT Male



Gas	Part No.
O ₂	6153
Air.....	6253
Vac*	6353

DISS Female



Gas	Part No.
Hand Tight O ₂	6158
Hand Tight Air.....	6258
Hand Tight Vac*	6358

*Vacuum primary check only

CUSTOM QUICK CONNECTS AND COUPLERS

OHMEDA FEMALE X OHMEDA FEMALE



Gas	Part #	Price
Oxygen.....	SP-105
Air.....	SP-105A
Vacuum	SP-105V

OHMEDA FEMALE COUPLER X DISS FEMALE HAND TIGHT WITH 90° ELBOW



Gas	Part #	Price
Oxygen.....	SP-51-OHF
Air.....	SP-51-AOHF
Vacuum	SP-51-VOHF

OHMEDA OFFSET "V" BLOCK



Gas	Part #	Price
Oxygen.....	SP-211
Air.....	SP-211A
Vacuum	SP-211V

CHEMETRON FEMALE X CHEMETRON FEMALE



Gas	Part #	Price
Oxygen.....	SP-106
Air.....	SP-106A
Vacuum	SP-106V

CHEMETRON FEMALE COUPLER X DISS FEMALE HAND TIGHT WITH 90° ELBOW



Gas	Part #	Price
Oxygen.....	SP-53CHF
Air.....	SP-53ACHF
Vacuum	SP-53VCHF

CHEMETRON OFFSET "V" BLOCK




Gas	Part #	Price
Oxygen.....	SP-210
Air.....	SP-210A
Vacuum	SP-210V

CHEMETRON FEMALE COUPLER X OHMEDA MALE QUICK CONNECT



Gas	Part #	Price
Oxygen.....	SP-129
Air.....	SP-129A
Vacuum	SP-129V

CHEMETRON MALE QUICK CONNECT X OHMEDA FEMALE COUPLER



Gas	Part #	Price
Oxygen.....	SP-130
Air.....	SP-130A
Vacuum	SP-130V

PURITAN BENNETT MALE X OHMEDA FEMALE COUPLER



Gas	Part #	Price
Oxygen.....	SP-209
Air.....	SP-209A
Vacuum	SP-209V

PURITAN BENNETT MALE X CHEMETRON FEMALE COUPLER



Gas	Part #	Price
Oxygen.....	SP-208
Air.....	SP-208A
Vacuum	SP-208V

OHMEDA VERTICAL "T" BLOCK



Gas	Part #	Price
Oxygen.....	SP-199
Air.....	SP-199A
Vacuum	SP-199V

CHEMETRON VERTICAL "T" BLOCK



Gas	Part #	Price
Oxygen.....	SP-196
Air.....	SP-196A
Vacuum	SP-196V

RC MEDICAL CAN ASSEMBLE MOST CUSTOM CONFIGURATIONS

PLEASE EMAIL: QUESTIONS@RCMEDICAL.COM WITH DRAWINGS OR REQUIREMENTS

DISS FITTINGS

DISS 90° Elbow Male x Female



Oxygen	Air	Vac
#0705	#SP-93	#SP-94

DISS Male x NPT Male



	Oxygen	Air	Vac
1/8"	#0701...	#0731...	#0751...
1/4"	#0703...	#0733...	#0753...

DISS Male x DISS Male



Oxygen	Air	Vac
#0704	#SP-91	#SP-92

DISS Male x 1/8" NPT Male with Check Valve



	Oxygen	Air	Vac
1/8"	#0712...	#0742...	#0762...
1/4"	#0713...	#0743...	#0763...

DISS Female x DISS Female



Oxygen	Air	Vac
#0104	#SP-88	#SP-89

DISS Male x DISS Female



Oxygen	Air	Vac
#SP-95	#SP-96	#SP-97

DISS Female x 1/8" NPT Male



Oxygen	Air	Vac
#0110	#0209	#0309

DISS Male x 1/8" NPT Female



Oxygen	Air	Vac
#0702	#SP-98	#SP-99





50 PSI
Auxiliary
Valve

OPTO
(OXYGEN)

APTO
(AIR)

Thread Comparisons Guide:



(Measurements given are approximate.)

For Oxequip, Medstar and Schrader
quick connects

PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM



NEB DRIVE RESCUE FLOW

An Economical adapter that delivers a preset flow from standard wall outlets.
It's a convenient, dependable way to administer aerosols.

DISS FEMALE



O.	6 to 7 LPM	8116
O.	7 to 8 LPM	8117
O.	8 to 9 LPM	8118
O.	10 to 11 LPM	8119
Air	6 to 7 LPM	8216
Air	7 to 8 LPM	8217
Air	8 to 9 LPM	8218
Air	10 to 11 LPM	8219
Rescue Flow	14 to 16 LPM	8115
O.	24 to 26 LPM	8114

OHMEDA



O.	6 to 7 LPM	8126
O.	7 to 8 LPM	8127
O.	8 to 9 LPM	8128
O.	10 to 11 LPM	8129
Air	6 to 7 LPM	8226
Air	7 to 8 LPM	8227
Air	8 to 9 LPM	8228
Air	10 to 11 LPM	8229
Rescue Flow	14 to 16 LPM	8125
O.	24 to 26 LPM	8124

CHEMETRON




O.	6 to 7 LPM	8136
O.	7 to 8 LPM	8137
O.	8 to 9 LPM	8138
O.	10 to 11 LPM	8139
Air	6 to 7 LPM	8236
Air	7 to 8 LPM	8237
Air	8 to 9 LPM	8238
Air	10 to 11 LPM	8239
Rescue Flow	14 to 16 LPM	8135
O.	24 to 26 LPM	8134

PURITAN-BENNETT



O.	6 to 7 LPM	8146
O.	7 to 8 LPM	8147
O.	8 to 9 LPM	8148
O.	10 to 11 LPM	8149
Air	6 to 7 LPM	8246
Air	7 to 8 LPM	8247
Air	8 to 9 LPM	8248
Air	10 to 11 LPM	8249
Rescue Flow	14 to 16 LPM	8145
O.	24 to 26 LPM	8144

To order Neb Drive or Rescue Flow with shut off valve, please add an "S" to the end of the part number.
The price is more.



Shut off valve allows you to leave neb drive in wall outlet.

REPLACEMENT RESPIRATORY OXYGEN SENSORS

Replaces: Ohmeda
6050-0004-110



Part # SS-20A

OEM Compatible with:
Analytical Ind.: PSR-11-915-4
C&S: G0-100
Ceramatec: CAG-10
City Technology: MOX-7
Envitec: OOM110
Datex-Ohmeda: 6050-0004-110
Integral Process: 65008
Maxtec: MAX-10, R112P10
Mercury Medical: 10-103-10
Nuova: E-15/4
Pacifitech: PT-20A
Sensidyne: SV-20A
Teledyne: R-30MED
Ventrex: V-20A
Vlamed: R-30V
VTI: 103105



Part # SS-03

OEM Compatible with:
Analytical Ind.: PSR-11-77
C&S: G0-070
Ceramatec: CAG-17
City Technology: MOX-4
Datex-Ohmeda: 6051-0000-219
Envitec: OOM105
Hamilton: HM-03
Hudson RCI: V-03, 5558
Integral Process: 65013
Maxtec: MAX-17, R116P10
Mercury Medical: 10-103-02
MSA: 804674
Newport: OM 100, VM 200
Nuova: E-77/O
Sensidyne: SV-03
SensorTech: ST-03
Teledyne: T-7
Ventrex: V-03
Vlamed: T-7
VTI: 103102

Fits ALL MiniOX!



Part # SS-11A

OEM Compatible with:
Air Shields: 6735142
Alphamed: R-17
Analytical Ind.: PSR-11-917-J
Smiths Medical > BCI: R-17
BioMed Devices: M-2, M-10
Bird Products: 66030
C&S: G0-030
Ceramatec: CAG-13
Datex-Ohmeda: 6051-0000-217
 Draeger: 6735142
Envitec: OOM103-1, OOM202-1
Hamilton: HM-11
Hudson RCI: 5567, V-11
Invivo: 9445
Lifecare: PLV-100, 37-079
Maxtec: MAX-13, R115P10
Mercury Medical: 10-103-05
MSA: 406931
Oxitron: 3310
Pacifitech: PT-11A
Medtronic > Puritan Bennett: 4-020933-00
Respironics: 37-079
Sensidyne: SV-11A
SensorTech: ST-11
Teledyne: R-17
Vlamed: R-41V
VTI: 100580



Part # SS-25A

OEM Compatible with:
Analytical Ind.: PSR-11-917-J1
C&S: G0-250
Maxtec: MAX-250ESF, R115P85
Mercury Medical: 10-103-08
Nuova: E-25/O
Pacifitech: PT-25A
Sensidyne: SV-25A
VTI: 103104



Part # SS-12A

OEM Compatible with:
Analytical Ind.: PSR-11-917M
C&S: G0-040
Ceramatec: CAG-12
City Technology: MOX-9
Datex-Ohmeda: 6051-0000-222
Envitec: OOM202
Hamilton Medical: 396008, 396009, HM-12
Hudson RCI: 5566, V-12
Integral Process: 65016
IMT Medical: PF300
Maxtec: MAX-12, R109P45
Mercury Medical: 10-103-06
MSA: 472062, 470644
Nuova: E-17/M
Pacifitech: PT-12A
Medtronic > Puritan Bennett: 4-072214-00
Siemens: 9004979-E347E, 6419332-E380E
Sensidyne: SV-12A
SensorTech: ST-01
Spacelabs: 90513
Taema: YR049700
Teledyne: R-24 MED, B-62321
Ventrex: V-12A
Vlamed: R-24Med
VTI Vascular Technology: 103100

ALL SENSORS

Unconditional 12 month Warranty!
Exceptional Quality and Performance!
Canned for increased shelf life and freshness!



FOR MORE SENSORS
PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM

MINIOX 200E OXYGEN ANALYZER




MiniOX 200E Oxygen Analyzer

Part # 600210.....
Comes complete with sensor, cable, tee and AA batteries
Replacement Sensor # 606940

The MiniOX 200E is best used anywhere an inline analyzer is needed, when alarms are not deemed important. It is especially useful in conjunction with older ventilators or oscillators. Designed with simplicity and portability in mind, this analyzer quickly adapts for use in tubing applications and is ideal for checking oxygen percent accuracy in line or on concentrations. Its patented ergonomic design fits comfortably in your hand and features one-touch on/off and calibration. The MiniOX 200E will calibrate to room air (20.9%) or 100% Oxygen. It is ideal for applications needing an external sensor.

- ① Large LCD – Easy-to-read characters
- ② Calibration Reminder – Reminds once a week
- ③ One-Touch Calibration – Auto calibrate to room air or 100% Oxygen
- ④ On/Off Button – One-touch on and off
- ⑤ Patented Ergonomic Design – Fits comfortably in your hand
- ⑥ External Sensor Cable – 10 ft. extended, 3 ft. retracted
- ⑦ Low Battery Indicator – 5,000+ Hours of battery life
- ⑧ Rubber Gasket – Prevents moisture from contact board
- ⑨ External Sensor – MiniOX 250E, 0-100% Oxygen, 24 month warranty



Approximately 5,000 hours with continuous use

MISCELLANEOUS

PLASTIC BAGS (DISPOSABLE)



PATIENT SET-UP (Cotton Drawstring)

Part No.	Description	Quantity	Price
0049	8 x 10	500	
0059	11 3/4 x 15	500	
XL-S2	16 x 18 x 3	250	

EQUIPMENT COVERS

Dust covers designed to protect ventilators and other equipment. Comes on a roll.

Part No.	Description	Quantity	Price
PLC-60	28 x 22 x 58 (clear)	55	
PLC-60/B	28 x 22 x 58 (buff)	55	
0084	28 x 22 x 56 (clear***)	50	
7083	28 x 22 x 56 (yellow)	50	
0126	48 x 25 x 42 (red)	50	

***printed, "CLEAN, READY FOR USE"

Main Flow Bacterial / Viral Filter



- 22mm ID x 15mm ID/22mm OD
- Bacterial filtration efficiency 99.999+%
- Viral filtration efficiency 99.99+%
 - Dead space: 42mL
- Flow resistance: approximately 1.5 cmH₂O at 60 LPM
- Individually packaged 50 to a case

Part Number 1605
per case of 50

ANESTHESIA ASSOCIATES PRESSURE MANOMETER



Part#	Description	Price
199	Gauge with 1/4" NPT	
274	Gauge with 5-55 cm H ₂ O adjustable pressure relief valve and "T" for supply tubing.	

Part#	Description	Price
199-MAX	Gauge with Pointer	
274-MAX	274 with Pointer	

Manometer has range of -40 to +80 cm H₂O in a 2 1/2" diameter chrome plated case with 1/2" NPT bottom connection. Has front zero correcting screw.



199-MAX

MAX is a resettable low friction pointer to capture maximum pressure or vacuum reading

POSEY CUFFLATOR™

ENDOTRACHEAL TUBE INFLATOR & MANOMETER



An accurate, convenient tracheal cuff inflator and manometer. No more syringes, electricity or complicated connections needed to check cuff pressure regularly. Air vent button and inflator bulb let you vary cuff pressure quickly and easily. Clip on back allows easy positioning. Gauge shows recommended pressure range in cm H₂O.

#8199 Cufflator with instrument case.....

FOR MORE PRODUCTS AND SPECIALS
PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM

#2200 Flexible Connector (Silicone Rubber)

22mm I.D. x 22mm I.D.

This adapter was designed almost identical to the ones listed below:

INTEC 2200 MARQUEST 172200
BENNETT 3443 AM SYSTEMS 1139
Instrumentation Ind. 101-00-3443



Part Number 2200
per case of 50
each

MISCELLANEOUS

Timeter Airdyne Hospital Grade Air Compressors



Uniquely designed for bed-side use. Special material absorbs most noise for quiet performance. Audible and visual high temperature and low pressure alarms with DC powered six-volt backup. Vibration-free gauge shows diminished pressure, impaired drying or exceeded flow capacity. Circuit breaker and thermal overload protector included. Continuous flow of 64 LPM at 50 psi.

Part Number TE20



Designed for bedside use. Quiet, vibration-free performance. Provides condensation-free compressed air to respiratory therapy equipment. Vibration-free gauge shows diminished pressure, impaired drying or exceeded flow capacity. Circuit breaker and thermal overload protector included. Continuous flow of 42 LPM at 50 psi.

Part Number TE35

DeVilbiss Pulmo-Aide® Medication Nebulizer



3 PRONG HOSPITAL GRADE PLUG

Durable AC-powered medication nebulizer with Pulmo-Neb® 60-day disposable nebulizer and tubing. Unit weighs 7 lbs. and has 9 - 16 LPM range. 5 year warranty. **115 vac**

DV4650D

BLENDERS - PRECISION MEDICAL

Air/O₂



Its unique. It is the only unit available that allows the clinician to turn the bleed control to ON or OFF providing quieter operation. This feature is especially useful in the NICU and PICU.

Part #	Flow	Rate	Price
5200	High Flow	15 - 120 lpm	
5300	Low Flow	3 - 30 lpm	

Heliox



Its unique. It is the only free-standing pneumatic HeliO, blending device on the market today. The ability to turn the bleed control collar ON or OFF provides ± 3% F102 accuracy as well as quieter operation and maximal Heliox tank durations.

Part #	Flow	%	Price
5570	High Flow	70 - 30	
5580	High Flow	80 - 20	
5470	Low Flow	70 - 30	
5480	Low Flow	80 - 20	

PRECISION MEDICAL EASY AIR COMPRESSOR



EasyAir Compressor is a high performance portable medical air compressor designed to reliably supply compressed air 24 hours a day. EasyAir is ideal for continuous use in any clinical environment. Its sound dampened case is easy to clean, allows quiet operation and is durable enough for years of heavy use.

PM15P Pressure Gauge
PM15F Flow Gauge

RC MEDICAL SELLS A WIDE VARIETY OF RESPIRATORY SUPPLIES AND EQUIPMENT. PLEASE VISIT OUR WEBSITE AT: WWW.RCMEDICAL.COM

REGULATOR AND FLOWMETER ACCESSORIES




Nut and Nipple (Christmas Tree)

Part #	Description	Price
7-4A*	One piece Green	cs/50
7-5A*	Two piece Green	cs/50
7-5C*	Two piece Yellow	cs/50
7-5H*	Two piece Clear	cs/50

Plastic - For single Patient Use (Above)*

0159	Metal Nut/Nipple	each
0159C	Metal w/ 6" chain	each

CYLINDER WRENCHES FOR SMALL AND LARGE CYLINDERS



Part #	Description	Price
7-92A	Plastic Small Wrench, BLACK	cs/10
7-92B	Plastic Small Wrench, GREEN	cs/10
7-92D	Plastic Small Wrench, PINK	cs/10
66071	Plastic Small Hand Wheel	cs/10
MCW-2B	Metal Small Cylinder Wrench	cs/10
MCW-2BC	Metal Small Wrench with Chain	cs/10
5080	Hudson Small Metal Wrench	each
WR-3	Metal Wrench Small & Large	each



REPLACEMENT GAUGES FOR OXYGEN REGULATORS

2" Bottom Port
1/4" NPT

2510 2530

2543

1 1/2" Back Port
1/8" NPT

2501

2521

Low Profile
3000 PSI

MJ-C1

Part #	Description	Price
2510	15 lpm - 2" Bottom Port	
2530	Contents (0-4000 PSI) 2" Bottom	
2543	100 PSI 2" Bottom Port	
2501	15 lpm - 1 1/2" Back Port	
2521	Contents (0-4000 PSI) 1 1/2" Back	
2500	15 lpm - 1 1/2" Bottom Port	
2520	Contents (0-4000 PSI) 1 1/2" Bottom	
MJ-C1	Low Profile Contents (3000 psi)	
AIR, CO2 GAUGES		
MG-C11	AIR, CO2 (0 - 100 PSI) 2" Bottom	
MG-C9	AIR, CO2 (0-4000 PSI) 2" Bottom	
MG-C13	AIR, CO2 (0-4000 PSI) 1 1/2" Bottom	

Yoke Washer For Small Cylinder Regulators (CGA 870)



2544-2

Part #	Description	Price
2544-2	aluminum w/ rubber center	each
2544-2C	aluminum w/ rubber center	cs/50

LITERMETER

Every therapist should have one for checking the accuracy of regulators and flowmeters.



Part #	LPM	Price
180-5004	0-2.5	
180-5005	0-8	
180-5006	2-15	

ALL IN STOCK ORDERS PLACED BY 2 PM (EST) SHIPPED THE SAME DAY.

"Y" BLOCKS

Ohmeda Style

(Ohio)



2152	O ₂	
2252	Air	
2352	Vac	

Chemetron Style

(NCG)



3152	O ₂	
3252	Air	
3352	Vac	

REUSABLE NUT AND NIPPLE ADAPTER (CHRISTMAS TREE)

**HEX NUT
TUBING NIPPLE**



0159C

0159

0159.....

0159C.....

**HAND TIGHT
TUBING NIPPLE**



0154C

0154

0154.....

0154C.....



Main Flow Bacterial / Viral Filter

- 22mm ID x 15mm ID/22mm OD
- Bacterial filtration efficiency 99.999+%
- Viral filtration efficiency 99.99+%
 - Dead space: 42mL
- Flow resistance: approximately 1.5 cmH₂O at 60 LPM
- Individually packaged 50 to a case

Part Number 1605
 per case of 50

**YOKE WASHER
 (FOR SMALL CYLINDER
 REGULATORS)**



2544-2

2544-2C /cs. 50

**ADJUSTABLE HEIGHT
 SMALL CYLINDER CART
 FOR D & E CYLINDERS**



#PX-1060.....