

RESPIRATORY CATALOG

EVERYDAY AND SPECIALTY PRODUCTS



Order online at: www.rcmedical.com



RC MEDICAL, INC.

Phone: 888-431-2727

Fax: 888-337-0443

P.O. Box 833 • 284 Merrow Road • Tolland, CT 06084

HEALTHCARE FACILITIES HAVE AN OPEN 30 DAY ACCOUNT

NO MINIMUM ORDER - PURCHASE ORDERS ACCEPTED

F.O.B. Origin. • Prices subject to change without notice.
 30% restocking charge prior to 30 days. • Please call for a RGA Number.
 Products may NOT be returned after 30 days.

CURRENT PRICING CAN BE FOUND ON WEBSITE

INDEX

ADAPTERS

Disposable 19
 DISS 27
 Female Couplers..... 25
 Male Quick Connects 24
 NPT Fittings 22
 Y Blocks and T Blocks..... 23
 BAGS, EQUIPMENT COVERS 29
 COMPRESSORS 30
 CYLINDER CARTS AND STANDS 9 - 11
 FITTINGS (under ADAPTERS)
 FLOWMETERS 5 - 8
 FLOWMETER ACCESSORIES 31

HOSES 17, 18
 MISCELLANEOUS 29, 30
 OXYGEN ANALYZERS AND SENSORS..... 28
 RAIL SYSTEMS..... 16
 REGULATORS
 Air 15
 CO2 15
 Heliox 15
 Oxygen 12, 13
 Oxygen - Specialty..... 14
 Suction..... 2 - 4
 RESPIROMETERS 21

Amvex
 Continuous / Intermittent
 Suction Regulators



SPECIAL PRICING

WHILE QUANTITIES

LAST!



Ohmeda
Wall Adapter



Part #
LQ-21

Chemetron
Wall Adapter



Part #
LQ-22

Puritan
Wall Adapter





Part #
LQ-23

DISS
Wall Adapter




Part #
LQ-24

 All suction regulators above come standard with tubing nipple outlet.

 For DISS Male outlet, add a "D" to end of part number.

Toll Free Phone
 888-431-2727
 1-860-871-4156

Toll Free Fax
 888-337-0443
 1-860-870-0443


**SUCTION REGULATORS – by Precision Medical
CONTINUOUS AND CONTINUOUS/INTERMITTENT**

**MADE IN USA
10 YEAR WARRANTY**

**USER FRIENDLY LARGE KNOB
EASY TO READ FACE**



2 Mode Continuous

Wall Adapter
Inlet

Tubing Nipple
Outlet

DISS Male
Outlet

DISS Female
Ohmeda
Chemetron
Puritan



PM3003
PM3004
PM3005
PM3008



PM3014
PM3015
PM3016
PM3019



3 Mode Continuous/Full line

Wall Adapter
Inlet

Tubing Nipple
Outlet

DISS Male
Outlet

DISS Female
Ohmeda
Chemetron
Puritan



PM3103
PM3104
PM3105
PM3108



PM3114
PM3115
PM3116
PM3119



3 Mode Continuous/Intermittent

Wall Adapter
Inlet

Tubing Nipple
Outlet

DISS Male
Outlet

DISS Female
Ohmeda
Chemetron
Puritan



PM3303
PM3304
PM3305
PM3308



PM3314
PM3315
PM3316
PM3319

WALL ADAPTERS (INLET)



DISS Female
w/hand nut



Ohio (Ohmeda)




Chemetron (NCG)



Puritan

VACUTRON SUCTION REGULATORS by Allied

The new, sleekly contoured Vacutron® is a precision engineered instrument that combines the most desirable features for suction regulators. Smaller in size, but big on benefits, it offers advantages that include:




Continuous Suction Regulator

Part	#Inlet	Outlet
22-13-1106	Ohmeda	DISS Male
22-13-1107	DISS Female	DISS Male
22-13-1108	Chemetron (NCG)	DISS Male
22-13-1105	Puritan-Bennett	DISS Male
22-13-1206	Ohmeda	Tubing Nipple
22-13-1207	DISS Female	Tubing Nipple
22-13-1208	Chemetron (NCG)	Tubing Nipple
22-13-1205	Puritan-Bennett	Tubing Nipple

Continuous/Intermittent Suction Regulator

22-15-1106	Ohmeda	DISS Male
22-15-1107	DISS Female	DISS Male
22-15-1108	Chemetron (NCG)	DISS Male
22-15-1105	Puritan-Bennett	DISS Male
22-15-1206	Ohmeda	Tubing Nipple
22-15-1207	DISS Female	Tubing Nipple
22-15-1208	Chemetron (NCG)	Tubing Nipple
22-15-1205	Puritan-Bennett	Tubing Nipple

Collection Canisters



Disposable Collection Canisters

Size	Outlet	Units/Case	Part Number
2400ML	DISS	36/CASE	20-08-0001
2400ML	STEM	36/CASE	20-08-0002
1500ML	STEM	48/CASE	20-08-0003
1500ML	DISS	48/CASE	20-08-0004
2400ML	DISS	36/CASE w/tubing	20-08-0005
2400ML	STEM	36/CASE w/tubing	20-08-0006
1500ML	STEM	42/CASE w/tubing	20-08-0007
1500ML	DISS	42/CASE w/tubing	20-08-0008

Suction Canister Holder and Wall Bracket



White Plastic Holder
20-02-0150



Wall Bracket
22-11-0006

For suction accessories please visit our website at:

www.rcmedical.com

Tubing
Fittings
Vac Trap

Stands
Bacteria Filters
Wye Blocks

PRECISION MEDICAL CHROME BODY FLOWMETER

**5 Year
Warranty**



4MFA1001



1MFA9001



3MFA1001



0 - 1 LPM PEDIATRIC OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER
4MFA1001	None - 1/8" NPT Female
4MFA1003	DISS FEMALE w/Hand Tight
4MFA1006	CHEMETRON (NCG)
4MFA1005	OHIO/OHMEDA
4MFA1008	PURITAN BENNETT

0 - 3 LPM PEDIATRIC OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER
1MFA3001	None - 1/8" NPT Female
1MFA3003	DISS FEMALE w/Hand Tight
1MFA3006	CHEMETRON (NCG)
1MFA3005	OHIO/OHMEDA
1MFA3008	PURITAN BENNETT

0 - 15 LPM OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER
1MFA1001	None - 1/8" NPT Female
1MFA1003	DISS FEMALE w/Hand Tight
1MFA1006	CHEMETRON (NCG)
1MFA1005	OHIO/OHMEDA
1MFA1008	PURITAN BENNETT

0 - 15 LPM OXYGEN FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER

PART NUMBER	ADAPTER
1MFA1003PTO	DISS FEMALE w/Hand Tight
1MFA1006PTO	CHEMETRON (NCG)
1MFA1005PTO	OHIO/OHMEDA
1MFA1008PTO	PURITAN BENNETT
1MFA100MPTO	MEDSTAR

0 - 70 LPM HI FLOW OXYGEN FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER
1MFA8001	None - 1/8" NPT Female
1MFA8003	DISS FEMALE w/Hand Tight
1MFA8006	CHEMETRON (NCG)
1MFA8005	OHIO/OHMEDA
1MFA8008	PURITAN BENNETT

0 - 70 LPM HI FLOW AIR FLOWMETER WITH ADAPTER

PART NUMBER	ADAPTER
1MFA9001	None - 1/8" NPT Female
1MFA9003	DISS FEMALE w/Hand Tight
1MFA9006	CHEMETRON (NCG)
1MFA9005	OHIO/OHMEDA
1MFA9008	PURITAN BENNETT

SELECT FLOWMETER

15 LPM OXYGEN FLOWMETER WITH ADAPTER WITH FLOW SELECTOR BUILT IN

PART NUMBER	ADAPTER
3MFA1001	None - 1/8" NPT Female
3MFA1003	DISS FEMALE w/Hand Tight
3MFA1006	CHEMETRON (NCG)
3MFA1005	OHIO/OHMEDA
3MFA1008	PURITAN BENNETT

FOR MORE FLOWMETERS, PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM


FLOW SELECTOR



Flow Selector Features

- Give nebulizer treatments without disconnecting the humidifier
- Reduces the changeover time to seconds
- No more loose tubing laying around
- Eliminates tubing nipple adaptors (It has 2 of its own.)
- Fits all flowmeters that have DISS male outlet fittings
- Lightweight

Part #..... 1000



INCREMENT SCALES

0-200cc, 0-1LPM, 0-3 LPM,
0-70 LPM, 0-15 LPM

PRECISION MEDICAL CHROME BODY 15 LPM FLOWMETERS

OXYGEN

Adapters and quick connects may be purchased separate on pages 24 and 25



1MFA1005

15 LPM OXYGEN FLOWMETER WITH ADAPTER

Part #	Adapter
1MFA1001	None - 1/8" NPT Female
1MFA1002	Diss Female
1MFA1003	Diss Female w/hand nut
1MFA1006	Chemetron (NCG)
1MFA1005	Ohio (Ohmeda)
1MFA1008	Puritan
1MFA100M	Medstar
1MFA1007	Oxequip



Diss Female



1MFA1006PTO

15 LPM OXYGEN FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER

Part #	Adapter
1MFA1001 PT0	None - 1/8" NPT Female
1MFA1002 PT0	Diss Female
1MFA1003 PT0	Diss Female w/hand nut
1MFA1006 PT0	Chemetron (NCG)
1MFA1005 PT0	Ohio (Ohmeda)
1MFA1008 PT0	Puritan
1MFA100M PT0	Medstar
1MFA1007 PT0	Oxequip



Diss Female with hand nut



Chemetron (NCG)



Y1MFA1006

15 LPM DUAL OXYGEN FLOWMETER WITH ADAPTER

Part #	Adapter
Y1MFA1001	None - 1/8" NPT Female
Y1MFA1002	Diss Female
Y1MFA1003	Diss Female w/hand nut
Y1MFA1006	Chemetron (NCG)
Y1MFA1005	Ohio (Ohmeda)
Y1MFA1008	Puritan
Y1MFA100M	Medstar
Y1MFA1007	Oxequip



Ohio (Ohmeda)



Puritan



Y1MFA1005PT2

15 LPM DUAL OXYGEN FLOWMETER WITH (2) 50 PSI AUXILIARY VALVES AND ADAPTER

Part #	Adapter
Y1MFA1001 PT2	None - 1/8" NPT Female
Y1MFA1002 PT2	Diss Female
Y1MFA1003 PT2	Diss Female w/hand nut
Y1MFA1006 PT2	Chemetron (NCG)
Y1MFA1005 PT2	Ohio (Ohmeda)
Y1MFA1008 PT2	Puritan
Y1MFA100M PT2	Medstar
Y1MFA1007 PT2	Oxequip



Medstar



Oxequip

CHAIN-RING ASSEMBLY



The Chain-Ring Assembly is used to permanently attach Flowmeters to the wall outlet.

Part #	Description
8004	6" Chain w/Closed Ring
8005	6" Chain w/Open Ring



50 PSI Auxiliary Valve

PRECISION MEDICAL CHROME BODY 15 LPM FLOWMETERS

AIR



1MFA2003



1MFA2008PTO



Y1MFA2003



Y1MFA2005PT2

15 LPM AIR FLOWMETER WITH ADAPTER

Part #	Adapter
1MFA2001	None - 1/8" NPT Female
1MFA2002	Diss Female
1MFA2003	Diss Female w/hand nut
1MFA2006	Chemetron (NCG)
1MFA2005	Ohio (Ohmeda)
1MFA2008	Puritan
1MFA200M	Medstar
1MFA2007	Oxequip

15 LPM AIR FLOWMETER WITH 50 PSI AUXILIARY VALVE AND ADAPTER

Part #	Adapter
1MFA2001 PT0	None - 1/8" NPT Female
1MFA2002 PT0	Diss Female
1MFA2003 PT0	Diss Female w/hand nut
1MFA2006 PT0	Chemetron (NCG)
1MFA2005 PT0	Ohio (Ohmeda)
1MFA2008 PT0	Puritan
1MFA200M PT0	Medstar
1MFA2007 PT0	Oxequip

15 LPM DUAL AIR FLOWMETER WITH ADAPTER

Part #	Adapter
Y1MFA2001	None - 1/8" NPT Female
Y1MFA2002	Diss Female
Y1MFA2003	Diss Female w/hand nut
Y1MFA2006	Chemetron (NCG)
Y1MFA2005	Ohio (Ohmeda)
Y1MFA2008	Puritan
Y1MFA200M	Medstar
Y1MFA2007	Oxequip

15 LPM DUAL AIR FLOWMETER WITH (2) 50 PSI AUXILIARY VALVES AND ADAPTER

Part #	Adapter
Y1MFA2001 PT2	None - 1/8" NPT Female
Y1MFA2002 PT2	Diss Female
Y1MFA2003 PT2	Diss Female w/hand nut
Y1MFA2006 PT2	Chemetron (NCG)
Y1MFA2005 PT2	Ohio (Ohmeda)
Y1MFA2008 PT2	Puritan
Y1MFA200M PT2	Medstar
Y1MFA2007 PT2	Oxequip

Adapters and quick connects may be purchased separate on pages 24 and 25



Diss Female



Diss Female with hand nut



Chemetron (NCG)



Ohio (Ohmeda)



Puritan



Medstar




Oxequip



50 PSI Auxiliary Valve

FIXED FLOW OXYGEN FLOWMETERS (1/2 - 15 LPM)

Position Increments (LPM)
Off, .50, 1, 2, 3, 4, 5, 6, 8, 10, 15



FOR MORE FLOWMETERS, PLEASE VISIT
OUR WEBSITE AT:
WWW.RCMEDICAL.COM







It's Durable - The lack of a flow tube and the durable metal body prevents breakage problems long associated with standard flowmeters.

It's Compact - The compact size allows the unit to fit into tight headwall areas and makes it virtually damage proof when other equipment is installed.

It's Accurate - The patented multi-spring technology insures accuracy at every flow setting without the need for a visible float as in Thorpe tube flowmeters

It's User-Friendly - The knurled flow adjustment knob and the large, easy-to-read numbers make flow adjustment easy.

It's Versatile - The unit is available for headwall outlets as well as for ceiling mount drops. Both air and oxygen models are available.

	Adapters	Side Connection	Top Connection
	1/8" NPT FE	7MFA1001	7MFA1001TC
	DISS Female	7MFA1003	7MFA1003TC
	Ohmeda	7MFA1005	7MFA1005TC
	Chemetron	7MFA1006	7MFA1006TC
	Puritan-Bennett	7MFA1008	7MFA1008TC



NEB BLOCK w/FLOWMETER

Features:

- The Neb Block allows a gas outlet to simultaneously or individually provide various respiration treatments without disconnecting any tubing or equipment.
- Provides easy to use On/Off switch to quickly access oxygen flow to a hand held nebulizer or resuscitation bag.
- Reliable fixed flow orifices provide six ranges of flow for aerosol therapy or resuscitation.



6-7 LPM 8-9 LPM 14-16 LPM
7-8 LPM 10-11 LPM

 Ohmeda	6-7 LPM	8MFA6065
	7-8 LPM	8MFA6075
	8-9 LPM	8MFA6085
	10-11 LPM	8MFA6105
	14-16 LPM	8MFA6155
 Chemetron (NCG)	6-7 LPM	8MFA6066
	7-8 LPM	8MFA6076
	8-9 LPM	8MFA6086
	10-11 LPM	8MFA6106
	14-16 LPM	8MFA6156
 DISS Female Hand Tight	6-7 LPM	8MFA6063
	7-8 LPM	8MFA6073
	8-9 LPM	8MFA6083
	10-11 LPM	8MFA6103
	14-16 LPM	8MFA6153
 Puritan Bennett	6-7 LPM	8MFA6068
	7-8 LPM	8MFA6078
	8-9 LPM	8MFA6088
	10-11 LPM	8MFA6108
	14-16 LPM	8MFA6158

SPECIALTY OXYGEN FLOWMETERS



15 LPM O2 Flowmeter for Chemetron Wall Outlet with 50 PSI-Chemetron Female Quick Connect
Part #..... SFO-5



15 LPM O2 Flowmeter for Chemetron Wall Outlet with 50 PSI-Chemetron Female Quick Connect
Part #..... SFO-2



Dual 15 LPM Oxygen Flowmeter with Chemetron Quick Connect Inlet and (2) Chemetron Couplers (Auxiliary Valves)
Part #..... SFO-12



15 LPM O2 Flowmeter for Ohmeda Wall Outlet with 50 PSI-Ohmeda Female Quick Connect
Part #..... SFO-4




15 LPM O2 Flowmeter for Ohmeda Wall Outlet with 50 PSI-Ohmeda Female Quick Connect
Part #..... SFO-1



Dual 15 LPM Oxygen Flowmeter with Ohmeda Quick Connect Inlet and (2) Ohmeda Couplers (Auxiliary Valves)
Part #..... SFO-11

SMALL CYLINDER CARTS (D&E)

Part #.....6105



Part #.....PX-1060



Part #.....6103-4

41" high
3" casters swivel w/brake



These four-wheeled units offer better control and easy maneuverability in and out of tight spaces. They require no tilting or balancing, giving you total control over the cylinder load.

Part # 1060-F



Holds 4

Part # 6205



Holds 2

LARGE CYLINDER CARTS AND STANDS (H,K)

Dual
6214



Single
6114



All models have 4 wheels.

2 Wheels Only



Part # HCC2-M



Part #.....610FW-E



Fiberglass Cylinder Bases
Part #.....HS-2

Part #

Description

6114

Single (Most Popular)

6114-PN

Single (Pneumatic Wheels)

6214

Dual

6214-PN

Dual (Pneumatic Wheels)

SMALL CYLINDER CARTS (D&E)

2-Wheel Cart



Part #	Cylinders Held	Dimensions (Height-Depth-Width)
6041	4	42" - 10" - 18"
6061	6	42" - 15" - 18"
6121	12	42" - 20" - 22"

With Sound Deadening Kit (below)

6041-SDK	4	42" - 10" - 18"
6061-SDK	6	42" - 15" - 18"
6121-SDK	12	42" - 20" - 22"

4-Wheel Cart



Part #	Cylinders Held	Dimensions (Height-Depth-Width)
6064	6	42" - 15" - 18"
6124	12	42" - 20" - 22"
6184*	18	42" - 30" - 16"
6244*	24	42" - 30" - 20"

With Sound Deadening Kit (below)

6124-SDK	12	42" - 20" - 22"
6184-SDK*	18	42" - 30" - 16"
6244-SDK*	24	42" - 30" - 20"

* -4" wheels

Multiple D & E Size Cylinder Stands With Wheels



Part #	Cylinders Held	Dimensions (Height-Depth-Width)
6040-PC	4	23" - 13" - 10"
6060-PC	6	23" - 12" - 15"
6120-PC	12	23" - 15" - 22"

With Sound Deadening Kit (below)

6040-PC-SDK	4	23" - 13" - 11"
6060-PC-SDK	6	23" - 12" - 15"
6120-PC-SDK	12	23" - 15" - 22"

ALL CARTS ABOVE ARE PAINTED GREEN

RED CARTS - ADD A "R" TO END OF PART NUMBER

YELLOW CARTS - ADD A "Y" TO END OF PART NUMBER

SOUND DEADENING KITS (for Carts Above)



Plastic Bumper Guards

Sound Deadening Matting

Anthony has developed an exclusive sound deadening system for their multiple cylinder carts, used to distribute cylinders in and around hospitals, rest homes, convalescent homes, etc. or where QUIET is essential.

These sound deadening kits include heavy-duty sound deadening matting and plastic bumper guards. This formula eliminates the rattling & clanging of steel & aluminum cylinders while in transport.

SOUND DEADENING KITS ALSO SOLD SEPARATELY
PLEASE VISIT OUR WEBSITE: WWW.RCMEDICAL.COM



10" x 14" Plastic Sign
Part # SP-51



10" x 14" Plastic Sign
Part # SP-52

SMALL CYLINDER STANDS AND HOLDERS (D&E)

Multiple D & E Size Cylinder Stands



Part #	Cylinders Held	Dimensions (Height-Depth-Width)
6040	4	20" - 13" - 11"
6060	6	20" - 12" - 15"
6120	12	20" - 15" - 22"
6180	18	20" - 15" - 32"
6240	24	20" - 20" - 32"
6400	40	20" - 40" - 25"

With Sound Deadening Kit (below)

6040-SDK	4	20" - 13" - 11"
6060-SDK	6	20" - 12" - 15"
6120-SDK	12	20" - 15" - 22"
6180-SDK	18	20" - 15" - 32"
6240-SDK	24	20" - 20" - 32"
6400-SDK	40	20" - 40" - 25"

Can be bolted to floor or free standing.

ALL STANDS ABOVE ARE PAINTED GREEN
 RED STANDS - ADD A "R" TO END OF PART NUMBER
 YELLOW STANDS - ADD A "Y" TO END OF PART NUMBER



MULTIPLE CYLINDER STANDS



CR-EPC6 6 CYLINDERS D OR E RACK 16" x 11" x 10"
 CR-EPC12 12 CYLINDERS D OR E RACK 22" x 18" x 10"

Bolts to floor or free standing.



For any solid surfaces,
 crash carts, more...
 Part # 610WM




For wheelchairs
 Part # 610WC

FOR MORE INFORMATION ON SMALL CYLINDER BASES,
 CYLINDER HOLDERS FOR BEDS, IV POLES, M6 CYLINDER CARTS AND STANDS,
 AND MORE, PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM

OXYGEN REGULATORS - LITERS PER MINUTE

Responsive Respiratory Fixed Flow Oxygen Regulator CGA 870 for Small Cylinder (D & E)




Part #120-3040
for small cylinder (CGA 870)

ALL BRASS

0, .5, 1, 1.5, 2, 3, 4, 5, 6, 8, 10, 15 settings



Part #21017
for small cylinder (870)



Part #120-1061C
for large cylinder (CGA 540)

0, 1, 2, 4, 6, 8, 10, 12, 15, 18, 20, 25 settings



Two High Pressure Outlets Also

Part #21025
1 - 25 LPM

Western Medica Oxygen Regulators

(Diaphragm Type)

6 Year Warranty



15 LPM

Single stage regulator for large cylinders. 15 liters per minute.
ORDER #M1-540-15 FG

With Gauge-Savers
ORDER #M1-540-15 FGP



M1-540-15 FM

Single stage flowmeter/regulator for large cylinders. 15 liters per minute.

ORDER #M1-540-15 FM

Single stage flowmeter/regulator for small cylinders. 15 liters per minute.

ORDER #M1-870-15 FM

MADE IN USA



15 LPM

Single stage regulator for small cylinders. 15 liters per minute.
ORDER #M1-870-15 FG

With Gauge-Savers
ORDER #M1-870-15 FGP



15 LPM

Single stage, swivel yoke regulator for small cylinders. 15 liters per minute.
ORDER #MR-870-15 FG

With Gauge-Savers
ORDER #MR-870-15 FGP

To order 8 or 5 liter regulators, add an 8 or a 5 to the end of our part number. The same price as our 15 liter regulators.

HIGH PRESSURE REGULATORS



Adjustable 0 - 100 psi with DISS Outlet

Gas	Part #	Cylinder Size
Oxy.	M1-540-PG	CGA 540 (Large)
Oxy.	M1-870-PG	CGA 870 (Small)
Air	M1-346-PG	CGA 346 (Large)
Air	M1-950-PG	CGA 950 (Small)



Preset 50 psi with DISS Outlet


Gas	Part #	Cylinder Size
Oxy.	M1A-540-P	CGA 540 (Large)
Air	M1A-346-P	CGA 346 (Large)
N2O	M1A-326-P	CGA 326 (Large)
CO2	M1A-320-P	CGA 320 (Large)



Preset 50 psi with DISS Outlet

Gas	Part #	Cylinder Size
Oxy.	M1A-870-P	CGA 870 (Small D&E)
Air	M1A-950-P	CGA 950 (Small D&E)
N2O	M1A-910-P	CGA 910 (Small D&E)
CO2	M1A-940-P	CGA 940 (Small D&E)

MADE IN USA



Preset 50 psi with Ohmeda Quick Connect

Gas	Part #	Cylinder Size
Oxy.	SRO-19	CGA 540 (Large)
Air	SRO-34	CGA 346 (Large)



Preset 50 psi with Chemetron Quick Connect

Gas	Part #	Cylinder Size
Oxy.	SRO-21	CGA 540 (Large)
Air	SRO-36	CGA 346 (Large)



Preset 50 psi with Ohmeda Quick Connect


Gas	Part #	Cylinder Size
Oxy.	SRO-30	CGA 870 (Small D&E)
Air	SRO-32	CGA 950 (Small D&E)



Preset 50 psi with Chemetron Quick Connect

Gas	Part #	Cylinder Size
Oxy.	SRO-29	CGA 870 (Small D&E)
Air	SRO-31	CGA 950 (Small D&E)

SPECIALTY OXYGEN REGULATORS



Fixed flow Regulator
With Ohmeda 50 psi
Outlet. .5 -25 lpm
0 - 25 LPM
Part # SRO-11



Fixed flow Regulator
With Chemetron 50 psi
Outlet. .5 - 25 lpm
Part # SRO-12



Fixed flow Regulator
With Puritan-Bennett
50 psi outlet. .5 - 25 lpm
Part # SRO-13



Oxygen Regulator/
Flowmeter 15 LPM,
Large Cylinder with
50 PSI Valve
Part # SRO-14



Oxygen Regulator with
Dual Chemetron Female
Couplers, 50 PSI, Large
Cylinder
Part # SRO-15



Oxygen Regulator with
Dual Ohmeda Female
Couplers, 50 PSI, Large
Cylinder
Part # SRO-16



Large Cylinder Oxygen
Regulator with (2)
15 LPM Flowmeters
Part # SRO-17



Large Cylinder Oxygen
Regulator with (2)
15 LPM Flowmeters and
(2) Auxiliary Valves
Part # SRO-18



Oxygen Regulator, Preset
50 PSI with Ohmeda
Quick Connect, Large
Cylinder
Part # SRO-19



Oxygen Regulator, Preset
50 PSI with Puritan
Bennett Quick Connect,
Large Cylinder
Part # SRO-20



Oxygen Regulator, Preset
50 PSI with Chemetron
Quick Connect, Large
Cylinder
Part # SRO-21



Oxygen Regulator with
15 LPM flowmeter and
50 PSI Ohmeda Quick
Connect, Large Cylinder
Part # SRO-22



Oxygen Regulator with
15 LPM Flowmeter and
50 PSI Chemetron Quick
Connect, Large Cylinder
Part # SRO-23



Oxygen Regulator, Preset
50 PSI with Chemetron
Quick Connect, Small
Cylinder (CGA 870)
Part # SRO-29



0 - 25 LPM Regulator
with Chemetron Female
50 PSI Outlet for Large
Cylinder (CGA 540)
Part # SRO-25



0 - 25 LPM Regulator
with Ohmeda Female
50 PSI Outlet for Large
Cylinder (CGA 540)
Part # SRO-26


Specialty AIR Regulators, Visit our Website at WWW.RCMEDICAL.COM

For Regulator Accessories See Page 31

HELIOX, AIR, CO₂ REGULATORS

Heliox Regulators

Part #	Ratio	For Cylinder size
2210-20	70 - 30	Small D and E CGA 890 Yoke
2211-20	80 - 20	Small D and E CGA 890 Yoke
2212-20	70 - 30	Large CGA 280 Nut/Nipple
2213-20	80 - 20	Large CGA 280 Nut/Nipple



0 - 20 LPM

WESTERN MEDICA MEDICAL REGULATORS (AIR or CARBON DIOXIDE)

Single-stage regulators preset to 50 psi or adjustable to 100 psi, or with a flowmeter. A sensitive neoprene diaphragm assembly provides optimum accuracy and sensitivity. Dual inlet filtering. Available for small or large cylinders.

AIR


CO₂

		AIR	CO ₂
Preset 50 PSI		M1A-346-P for large cylinder	M1-320-P for large cylinder
Adjustable 0 - 100 PSI		M1-346-PG for large cylinder	M1-320-PG for large cylinder
12 Liters Per Minute (lpm)		M1-346-15FM for large cylinder	M1-320-12FM for large cylinder
Preset 50 PSI		M1A-950-P for small cylinder	M1-940-P for small cylinder
Adjustable 0 - 100 PSI		M1-950-PG for small cylinder	M1-940-PG for small cylinder
12 Liters Per Minute (lpm)		M1-950-15FM for small cylinder	M1-940-12FM for small cylinder

EASTERN RAIL SYSTEMS

Simple Outlet Expansion allows your hospital to use its existing outlets, and to create the number and placement of outlets needed. Simply mount a length of our rail on the wall - attach the needed outlets to the rail with our thumb screw clamp - attach a length of high-pressure hose to the block, and run it back to your existing outlet.

- Easy external mounting - no need to shut down the unit.
- No more weight on outlets
- Shift the weight of humidifiers & suction bottles to the rail
- Eliminate overloading
- Reduce outlet leakage. Save time and money
- Increase efficiency
- Works with all piping systems
- Have outlets where you need them NOW!



KIT CONTAINS

Double Gas Outlet
Two Foot Rail
3 Foot Hose
Mounting Kit

CHEMETRON (Pictured)

Gas	Part #
Oxy	SP-241-O
Air	SP-241-A
Vac	SP-241-V
N2O	SP-241-N

OHMEDA

Gas	Part #
Oxy	SP-242-O
Air	SP-242-A
Vac	SP-242-V
N2O	SP-242-N

DISS

Gas	Part #
Oxy	SP-243-O
Air	SP-243-A
Vac	SP-243-V
N2O	SP-243-N

KIT CONTAINS

Two Single Gas Outlets
Four Foot Rail
Two 3 Foot Hoses
Mounting Kit, Wye Block

CHEMETRON (Pictured)

Oxy	SP-244-O
Air	SP-244-A
Vac	SP-244-V
N2O	SP-244-N

OHMEDA

Oxy	SP-245-O
Air	SP-245-A
Vac	SP-245-V
N2O	SP-245-N

DISS

Oxy	SP-246-O
Air	SP-246-A
Vac	SP-246-V
N2O	SP-246-N

KIT CONTAINS

Two Double Gas Outlets
Four Foot Rail
Two 3 Foot hoses
Mounting Kit, Wye Block

CHEMETRON (Pictured)

Oxy	SP-247-O
Air	SP-247-A
Vac	SP-247-V
N2O	SP-247-N



OHMEDA

Oxy	SP-248-O
Air	SP-248-A
Vac	SP-248-V
N2O	SP-248-N

DISS

Oxy	SP-249-O
Air	SP-249-A
Vac	SP-249-V
N2O	SP-249-N

**FOR MORE EASTERN RAIL PRODUCTS PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM**



CUSTOM HOSE ASSEMBLIES

Ordering is easy as 1, 2, 3 . . . 4


How to Order:

Western Medica manufactures custom hose assemblies that conform to CGS, and NFPA standards. All assemblies are pressure-tested. Select your requirement from the charts below.

Example: OXY-024-FNN-OHF

- 1) OXY Oxygen Service
- 2) 024 24" Long
- 3) FNN Inlet Connection
- 4) OHF Outlet Connection

Evacuation Hose (purple) price may be more, depending on your adaptor



Pricing is based on Inlet and Outlet Connection, Hose Length (Per Inch)

MINIMUM HOSE LENGTH — 1 FOOT

HOSE ORDERS FOR CUSTOM HOSES MAY NOT BE CANCELLED

Gas Service

OXY	Oxygen	Green
AIR	Medical Air	Yellow
SUC	Suction	White
N2O	Nitrous Oxide	Blue
EVC	Evacuation	Purple
CO ₂	Carbon Dioxide	Gray
NIT	Nitrogen	Black
GNR	General	Clear

Adapters pictured below may be ordered Separately.

When ordering adapters only, please specify gas.

Example: chm - oxy

All hose assemblies are 1/4" conductive (unless specified)

 #CHM Chemetron (NCG) Male	 #FNN Diss Female with Hex nut	 #FHT Diss Female with Hand Nut	 #OHM Ohmeda Male
 #MBA Diss Male	 #SHM Schrader Male	 #OXM Oxequip Male	 #PBM Puritan - Bennett Male
 #CHF Chemetron (NCG) Female	 #SHF Schrader Female	 #PBF Puritan - Bennett Female	 #OHF Ohmeda Female
 #DBA Demand Body Adapter	 #MP2 Male Pipe 1/8"	 #MP4 Male Pipe 1/4"	
 #FP2 Female Pipe 1/8"	 #FP4 Female Pipe 1/4"	 #BLK Blank (No Connection)	

REINFORCED MEDICAL HOSE



Hose only!

Conductive, Non-UL Listed
400 psi Working Pressure
1500 psi Burst Pressure

- MH4-04C..... 1/4" ID, Nitrous Oxide (Blue)
- MH4-08C..... 1/4" ID, Carbon Dioxide (Grey)
- MH4-12C..... 1/4" ID, Nitrogen (Black)
- MH4-16C..... 1/4" ID, Medical Air (Yellow)
- MH4-22C..... 1/4" ID, Suction (White)
- MH4-24C..... 1/4" ID, Oxygen (Green)
- MH5-22C..... 5/16" ID, Suction (White)

FERRULE



For 1/4" hose
.525 ID
#M-625-A

VICE-GRIPS

For 3/16", 1/4" and
5/16" ID Hose
C-5A



COMMON HOSE ASSEMBLIES

DISS FEMALE HANDTIGHT X DISS FEMALE HANDTIGHT



DESCRIPTION	PART #
5-Foot Oxygen	CH-100
5 Foot Air	CH-101
10 Foot Oxygen	CH-102
10 Foot Air	CH-103
15 Foot Oxygen	CH-104
15 Foot Air	CH-105
20 Foot Oxygen	CH-106
20 Foot Air	CH-107

DISS FEMALE HANDTIGHT X CHEMETRON MALE



DESCRIPTION	PART #
5-Foot Oxygen	CH-108
5 Foot Air	CH-109
10 Foot Oxygen	CH-110
10 Foot Air	CH-111
15 Foot Oxygen	CH-112
15 Foot Air	CH-113
20 Foot Oxygen	CH-114
20 Foot Air	CH-115

DISS FEMALE HANDTIGHT X CHEMETRON FEMALE



DESCRIPTION	PART #
5-Foot Oxygen	CH-116
5 Foot Air	CH-117
10 Foot Oxygen	CH-118
10 Foot Air	CH-119
15 Foot Oxygen	CH-120
15 Foot Air	CH-121
20 Foot Oxygen	CH-122
20 Foot Air	CH-123

CHEMETRON MALE X CHEMETRON FEMALE



DESCRIPTION	PART #
5-Foot Oxygen	CH-124
5 Foot Air	CH-125
10 Foot Oxygen	CH-126
10 Foot Air	CH-127
15 Foot Oxygen	CH-128
15 Foot Air	CH-129
20 Foot Oxygen	CH-130
20 Foot Air	CH-131

DISS FEMALE HANDTIGHT X PURITAN BENNETT MALE



DESCRIPTION	PART #
5-Foot Oxygen	CH-164
5 Foot Air	CH-165
10 Foot Oxygen	CH-166
10 Foot Air	CH-167
15 Foot Oxygen	CH-168
15 Foot Air	CH-169
20 Foot Oxygen	CH-170
20 Foot Air	CH-171

DISS FEMALE HANDTIGHT X OHMEDA MALE



DESCRIPTION	PART #
5-Foot Oxygen	CH-132
5 Foot Air	CH-133
10 Foot Oxygen	CH-134
10 Foot Air	CH-135
15 Foot Oxygen	CH-136
15 Foot Air	CH-137
20 Foot Oxygen	CH-138
20 Foot Air	CH-139

DISS FEMALE HANDTIGHT X OHMEDA FEMALE



DESCRIPTION	PART #
5-Foot Oxygen	CH-140
5 Foot Air	CH-141
10 Foot Oxygen	CH-142
10 Foot Air	CH-143
15 Foot Oxygen	CH-144
15 Foot Air	CH-145
20 Foot Oxygen	CH-146
20 Foot Air	CH-147

OHMEDA MALE X OHMEDA FEMALE



DESCRIPTION	PART #
5-Foot Oxygen	CH-148
5 Foot Air	CH-149
10 Foot Oxygen	CH-150
10 Foot Air	CH-151
15 Foot Oxygen	CH-152
15 Foot Air	CH-153
20 Foot Oxygen	CH-154
20 Foot Air	CH-155

PURITAN BENNETT MALE X PURITAN BENNETT FEMALE



DESCRIPTION	PART #
5-Foot Oxygen	CH-156
5 Foot Air	CH-157
10 Foot Oxygen	CH-158
10 Foot Air	CH-159
15 Foot Oxygen	CH-160
15 Foot Air	CH-161
20 Foot Oxygen	CH-162
20 Foot Air	CH-163

DISS FEMALE HANDTIGHT X PURITAN BENNETT FEMALE



DESCRIPTION	PART #
5-Foot Oxygen	CH-172
5 Foot Air	CH-173
10 Foot Oxygen	CH-174
10 Foot Air	CH-175
15 Foot Oxygen	CH-176
15 Foot Air	CH-177
20 Foot Oxygen	CH-178
20 Foot Air	CH-179

HUDSON RCI DISPOSABLE ADAPTERS AND CONNECTORS



No. 1667
Ventilator Monitoring
Adaptor Circuit

- Allows monitoring of ventilatory weaning parameters in intubated patients
- Includes oxygen nipple adaptor, aerosol tee, 6" CORR-A-FLEX tube, 15mm/22mm adaptor, m/m one-way valve and m/f one-way valve
- 50 per case



No. 1637
Thermometer
with Tee Adapter

- Combines No. 1635 Thermometer and No. 1077 Aerosol Tee Connector
- Packaged 50 to a case



No. 1743
Valved Neb Tee

- 22mm O.D. x 18mm I.D. x 22mm I.D. connections
- Tee may be left in-line when nebulizer is removed
- Port accepts most small volume nebulizer
- 30 per case



No. 1648
Aerosol Tee Connector

- 22mm O.D. x 15mm I.D./22mm O.D. x 22mm O.D.
- 50 per case



No. 1641
Elbow Adaptor

- Standard 22mm I.D. connector on vertical axis
- 22mm O.D. connection on horizontal axis
- Packaged 50 to a case



No. 1642
Pressure Line Adaptor

- 22mm I.D./O.D. connections
- Barbed elbow accepts 5 to 7mm I.D. tubing and swivels 360 degrees
- Packaged 50 to a case



No. 1623
Patient Elbow Adaptor

- 15mm O.D. x 15mm I.D./22mm O.D.
- 50 per case



No. 1644
One-Way Valve

- Standard 22 mm I.D./O.D. connections
- Valve offers low flow resistance and a capped monitoring port
- Packaged 50 to a case




No. 1638
Aerosol Tee Connector

- Horizontal ports of the tee are 22mm O.D.
- Vertical port is 18mm I.D./22mm O.D.
- Packaged 50 to a case



No. 1077
Aerosol Tee Connector

- Anti-spill design minimizes aspiration of liquids
- Horizontal axis has standard 22mm O.D. connections
- Vertical axis incorporates a standard 15mm I.D. connection and two 22mm O.D. concentric rings
- Packaged 50 to a case



No. 1639
Nebulizer Tee Connector

- Standard 22 mm I.D./O.D. connections
- Nebulizer port is 18mm I.D. which accepts most small volume nebulizers
- Packaged 50 to a case



No. 1421
Universal Cuff Adaptor

- Standard 22mm I.D. connections
- Packaged 50 to a case



No. 1422
Multi Adaptor

- Standard 15mm I.D./22mm O.D. connections
- Packaged 50 to a case



No. 1409
Large Tubing Connector

- Standard 22mm O.D. connections
- Packaged 50 to a case



No. 1640
Tubing Adaptor

- Standard 22mm I.D./O.D.
- Packaged 50 to a case



No. 1751
MDI Adapter

- 15mm O.D x 15mm I.D. connections for in-line placement in 15mm tubing
- 50 per case



No. 1659
MDI Adapter

- 22mm O.D x 22mm I.D. connections for in-line placement in 22mm tubing
- 50 per case



No. 1423
Tubing Adaptor

- 22mm I.D. x 5 - 7mm tapered nipple
- Packaged 50 to a case



No. 1669
Monitoring Adaptor

- 15mm O.D. x 15mm I.D./22mm O.D.
- Standard luer fitting
- 50 per case

DISPOSABLES

#2200 Flexible Connector (Silicone Rubber)

22mm I.D. x 22mm I.D.

This adapter was designed almost identical to the ones listed below:

INTEC 2200 MARQUEST 172200
 BENNETT 3443 AM SYSTEMS 1139
 Instrumentation Ind. 101-00-3443

Part Number 2200
Case of 50

Part Number 2200E
Eaches



Tapered Silicone Rubber Connector

27mm ID x 22mm ID

Part Number 1426
Case of 10



MDI Accessory by Hudson

The POCKET CHAMBER is equipped with a highly sensitive silicone valve that opens easily without restricting the air flow. This one-way valve helps keep patients from accidentally blowing medication out of the chamber, and may aid patients experiencing difficulty timing their inhalation.

Part Number 1001-10
Case of 10

Main Flow Bacterial / Viral Filter

- 22mm ID x 15mm ID/22mm OD
- Bacterial filtration efficiency 99.999+%
- Viral filtration efficiency 99.99+%
 - Dead space: 42mL
- Flow resistance: approximately 1.5 cmH₂O at 60 LPM
- Individually packaged 50 to a case

Part Number 1605
Case of 50



Supply Line Adapter

For extending smooth bore oxygen tubing
 Tapered for excellent gripping

Part Number 65100
Case of 200

Corrugated Tubing

- Continuous 100-foot roll of tubing with cuttable 6-inch sections
- Packaged in protective, convenient, plastic-lined dispenser box

Part Number 1680
roll / 100 FT



Ventilator Circuits

Adult Standard Ventilator Circuit, 72" long. Use on Puritan-Bennett 7200A; Servo 300, 900A, 900B, and 900C; Hamilton Veolar; Emerson; Engstrom Erica; PPG Irisa; Infracor Adult Star and Bird 8400ST

Part Number 1613
Case of 10



Voldyne Incentive Spirometer

Voldyne volumetric exerciser features Good - Better - Best flow window in a compact, ergonomic design with advanced, low work-of-breathing filter design and built in handle. 5000 ml

Part Number 8884719010
Case of 12



E-Z Wrap

Soft closed-cell foam tubes designed to prevent ear chafing. Light gray color Packed 50 pair per case.

Part Number 1016
Case of 50



Noseclip

Part Number 64019
LATEX FREE Package of 100



Aerosol Drainage Bags

Disposable Aerosol Drainage Bags with tee adaptor fits 22mm ID corrugated tubing. Bag has plastic chain for securing and convenient drainage port. 750 cc capacity. With wye adaptor.

Part Number 1742
Case of 50



RESPIROMETERS

MARK 8



1.4" Face

For measurement of minute volumes as well as tidal volumes. Two dials. Small dial records 0-1 liter. Large dial 0-100 liters.

Part Number 2750

MARK 14



1.4" Face

Specially suited for the measurement of tidal volumes and may also be used for measurement of minute volumes. Two dials. Small dial records 0-10 liters. Large dial records 0-1 liter.

Part Number 2751

HALOSCALE STANDARD



2.25" Face

The original design of Haloscale Respirometer incorporating concentric dials and an enlarged face for improved readability. Equally convenient for measuring both tidal and minute volumes. Large 57 mm dial diameter. Red/black hands and black scales on a white background with yellow highlight ring. Two concentric outer dials. First dial records 0-1 liter. Second dial records 0-100 liters.

Part Number 2768

REPAIRS

Includes cleaning, calibration and oiling. Parts are extra
All Models Repaired including:

**Wright • Ferraris • Haloscale • BOC • Fraser - Harlake
Ohmeda • Boehringer**

ELECTRONIC DIGITAL RESPIROMETER

A respirometer long overdue in this modern electronic age



ANESTHESIA ASSOCIATES

Part #	
295	Respirometer
Options	
296	6' Retractable Cable (extended length)
298	Carrying Case
299	Spare Flow Sensing Unit

MODEL 295 DIGITAL RESPIROMETER IS AN ELECTRONIC VERSION OF THE MECHANICAL RESPIROMETER. IT DOES THE SAME JOB PLUS

- Easier to read
- Saving on the initial purchase and maintenance
- More reliable
- Extremely durable
- No gears, shafts, or bearings to break down in the electronic unit
- Separates for ease of patient use, sterilization and cleaning

Model 295, a two piece detachable respirometer comes complete ready to use with 9v battery and one (1) 22 mm female x 22 mm female adapter.

**BROWSE OUR WEBSITE AT:
WWW.RCMEDICAL.COM**

Reconditioned Used Respirometers















Fully Serviced, Cleaned, Oiled and Calibrated

1 Year Warranty





Order #2799
(Case included)









NPT PIPE FITTINGS

<p>Female Tee</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7034 $\frac{1}{4}$" NPT..... 7035</p>	<p>Street Tee</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7032 $\frac{1}{4}$" NPT..... 7033</p>	<p>Cross</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7036 $\frac{1}{4}$" NPT..... 7037</p>	<p>45° Street Elbow</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7056 $\frac{1}{4}$" NPT..... 7057</p>
<p>Connectors</p>  <p>Part No. $\frac{1}{8}$" Female x $\frac{1}{8}$" Male 1017 $\frac{1}{8}$" Female x $\frac{1}{4}$" Male 7024 $\frac{1}{4}$" Female x $\frac{1}{8}$" Male..... 7025 $\frac{1}{4}$" Female x $\frac{1}{4}$" Male 7026</p>	<p>Hex Nipple</p>  <p>Part No. $\frac{1}{8}$" Male x $\frac{1}{8}$" Male..... 7020 $\frac{1}{4}$" Male x $\frac{1}{4}$" Male..... 7021 $\frac{1}{8}$" Male x $\frac{1}{4}$" Male..... 7022</p>	<p>Hose Barb Connector</p>  <p>Part No. $\frac{1}{8}$" Male x $\frac{1}{8}$" Dia. 7038 $\frac{1}{8}$" Male x $\frac{3}{16}$" Dia..... 7039 $\frac{1}{8}$" Male x $\frac{1}{4}$" Dia. 7040 $\frac{1}{4}$" Male x $\frac{1}{4}$" Dia. 7043</p>	<p>Coupler</p>  <p>Part No. $\frac{1}{8}$" Female x $\frac{1}{8}$" Female 7027 $\frac{1}{4}$" Female x $\frac{1}{4}$" Female 7028 $\frac{1}{4}$" Female x $\frac{1}{8}$" Female 7029</p>
<p>90° Elbow</p>  <p>Part No. $\frac{1}{8}$" Female x $\frac{1}{8}$" Female 7058 $\frac{1}{4}$" Female x $\frac{1}{4}$" Female 7059</p>	<p>Plug</p>  <p>Part No. $\frac{1}{4}$" NPT..... 7017 $\frac{1}{8}$" NPT..... 7019</p>	<p>Hose Barb Connector</p>  <p>Part No. $\frac{1}{8}$" Female x $\frac{1}{4}$" Dia. 7046 $\frac{1}{4}$" Female x $\frac{1}{4}$" Dia. 7049</p>	<p>Cap</p>  <p>Part No. $\frac{1}{4}$" NPT..... 7016 $\frac{1}{8}$" NPT..... 7018</p>
<p>Pipe Tee</p>  <p>Part No. $\frac{1}{8}$" NPT..... 7030 $\frac{1}{4}$" NPT..... 7031</p>	<p>Drop Ear</p>  <p>Part No. Tee $\frac{1}{8}$" F x $\frac{1}{8}$" F x $\frac{1}{8}$" F 7081 Elbow $\frac{1}{8}$" F x $\frac{1}{8}$" F..... 7082</p>	<p>Pipe Nipple</p>  <p>Part No. $\frac{1}{8}$" Male x $1\frac{1}{2}$" Long 7050 $\frac{1}{8}$" Male x 2" Long 7051 $\frac{1}{8}$" Male x $2\frac{1}{2}$" Long 7052 $\frac{1}{4}$" Male x $1\frac{1}{2}$" Long 7053 $\frac{1}{4}$" Male x 2" Long 7054 $\frac{1}{4}$" Male x $2\frac{1}{2}$" Long 7055 $\frac{1}{8}$" Male x $3\frac{1}{2}$" Long 7062</p>	<p>90° Elbow</p>  <p>Part No. $\frac{1}{8}$" F x $\frac{1}{8}$" M..... 7063 $\frac{1}{4}$" F x $\frac{1}{4}$" M..... 7065</p>

"Y" BLOCK ASSEMBLIES

Ohmeda Style (Ohio)	Chemetron Style (NCG)	Puritan-Bennett Style	Ohmeda Style "Y" Assembly—"Stubby"
			
2152..... O ₂ 2252..... Air 2352..... Vac	3152..... O ₂ 3252..... Air 3352..... Vac	4152..... O ₂ 4252..... Air 4352..... Vac	Compact Version 2162..... O ₂ 2262..... Air 2362..... Vac

T BLOCKS

OHMEDA Style			
	Gas Part # Oxy 2152-T Air 2252-T Vac 2352-T		Gas Part # Oxy 2152-TL Air 2252-TL Vac 2352-TL
	Gas Part # Oxy 2152-TR Air 2252-TR Vac 2352-TR		
CHEMETRON Style			
	Gas Part # Oxy 3152-T Air 3252-T Vac 3352-T		Gas Part # Oxy 3152-TL Air 3252-TL Vac 3352-TL
	Gas Part # Oxy 3152-TR Air 3252-TR Vac 3352-TR		

TRI-BLOCK ASSEMBLIES



Ohmeda Style
 SP-112 O₂
 SP-113 Air
 SP-114 Vac



Chemetron Style
 SP-115 O₂
 SP-116 Air
 SP-117 Vac

DISS "Y" ASSEMBLY WITH CHECK VALVES



(7211)

7211 O₂
 7212 Air
 7213 Vac

**FOR MORE "Y" BLOCK ASSEMBLIES
 PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM**

DISS "Y" WITH 2 OHMEDA COUPLERS



SP-214 O₂
 SP-214A Air
 SP-214V Vac

DISS "Y" WITH 2 CHEMETRON COUPLERS



SP-215 O₂
 SP-215A Air
 SP-215V Vac

CHANGE YOUR DISS WALL OUTLET TO QUICK CONNECTS

QUICK CONNECTS (FEMALE)

CHEMETRON STYLE CHECK UNITS (Female Couplers)

OHMEDA STYLE CHECK UNITS (Female Couplers)



1/8" NPT Female

DISS Male



1/8" NPT Female

DISS Male



Gas	Part No.	Gas	Part No.
O ₂	3150	O ₂	3155
Air.....	3250	Air.....	3255
Vac*.....	3350	Vac*.....	3355

Gas	Part No.	Gas	Part No.
O ₂	2150	O ₂	2155
Air.....	2250	Air.....	2255
Vac*.....	2350	Vac*.....	2355

1/4" NPT Male

DISS Female



1/4" NPT Male

DISS Female



Gas	Part No.	Gas	Part No.
O ₂	3153	Hand Tight O ₂	3158
Air.....	3253	Hand Tight Air.....	3258
Vac*.....	3353	Hand Tight Vac*	3358

Gas	Part No.	Gas	Part No.
O ₂	2153	Hand Tight O ₂	2158
Air.....	2253	Hand Tight Air.....	2258
Vac*.....	2353	Hand Tight Vac*	2358

PURITAN-BENNETT STYLE CHECK UNITS (Female Couplers)

SCHRADER STYLE CHECK UNITS (Female Couplers)

1/8" NPT Female

DISS Male



1/8" NPT Female

DISS Male



Gas	Part No.	Gas	Part No.
O ₂	4150	O ₂	4155
Air.....	4250	Air.....	4255
Vac*.....	4350	Vac*.....	4355

Gas	Part No.	Gas	Part No.
O ₂	6150	O ₂	6155
Air.....	6250	Air.....	6255
Vac*.....	6350	Vac*.....	6355

1/4" NPT Male

DISS Female



1/4" NPT Male

DISS Female



Gas	Part No.	Gas	Part No.
O ₂	4153	Hand Tight O ₂	4158
Air.....	4253	Hand Tight Air.....	4258
Vac*.....	4353	Hand Tight Vac*	4358

Gas	Part No.	Gas	Part No.
O ₂	6153	Hand Tight O ₂	6158
Air.....	6253	Hand Tight Air.....	6258
Vac*.....	6353	Hand Tight Vac*	6358

*Vacuum primary check only

CUSTOM QUICK CONNECTS AND COUPLERS

**OHMEDA FEMALE X
OHMEDA FEMALE**



Gas	Part #
Oxygen.....	SP-105
Air.....	SP-105A
Vacuum	SP-105V

**OHMEDA FEMALE COUPLER X DISS
FEMALE HAND TIGHT WITH 90° ELBOW**



Gas	Part #
Oxygen.....	SP-51-OHF
Air.....	SP-51-AOHF
Vacuum ...	SP-51-VOHF

OHMEDA OFFSET "V" BLOCK



Gas	Part #
Oxygen.....	SP-211
Air.....	SP-211A
Vacuum	SP-211V

**CHEMETRON FEMALE X
CHEMETRON FEMALE**



Gas	Part #
Oxygen.....	SP-106
Air.....	SP-106A
Vacuum	SP-106V

**CHEMETRON FEMALE COUPLER X DISS
FEMALE HAND TIGHT WITH 90° ELBOW**



Gas	Part #
Oxygen.....	SP-53CHF
Air.....	SP-53ACHF
Vacuum	SP-53VCHF

CHEMETRON OFFSET "V" BLOCK




Gas	Part #
Oxygen.....	SP-210
Air.....	SP-210A
Vacuum	SP-210V

**CHEMETRON FEMALE COUPLER X
OHMEDA MALE QUICK CONNECT**



Gas	Part #
Oxygen.....	SP-129
Air.....	SP-129A
Vacuum	SP-129V

**CHEMETRON MALE QUICK CONNECT
X OHMEDA FEMALE COUPLER**



Gas	Part #
Oxygen.....	SP-130
Air.....	SP-130A
Vacuum	SP-130V

**PURITAN BENNETT MALE X OHMEDA
FEMALE COUPLER**



Gas	Part #
Oxygen.....	SP-209
Air.....	SP-209A
Vacuum	SP-209V

**PURITAN BENNETT MALE X
CHEMETRON FEMALE COUPLER**



Gas	Part #
Oxygen.....	SP-208
Air.....	SP-208A
Vacuum	SP-208V

OHMEDA VERTICAL "T" BLOCK



Gas	Part #
Oxygen.....	SP-199
Air.....	SP-199A
Vacuum	SP-199V

CHEMETRON VERTICAL "T" BLOCK



Gas	Part #
Oxygen.....	SP-196
Air.....	SP-196A
Vacuum	SP-196V

RC MEDICAL CAN ASSEMBLE MOST CUSTOM CONFIGURATIONS

PLEASE EMAIL: QUESTIONS@RCMEDICAL.COM WITH DRAWINGS OR REQUIREMENTS

DISS FITTINGS

DISS 90° Elbow Male x Female



Oxygen	Air	Vac
#0705	#SP-93	#SP-94

DISS Male x NPT Male



	Oxygen	Air	Vac
1/8"	#0701	#0731	#0751
1/4"	#0703	#0733	#0753

DISS Male x DISS Male



Oxygen	Air	Vac
#0704	#SP-91	#SP-92

DISS Male x 1/8" NPT Male with Check Valve



	Oxygen	Air	Vac
1/8"	#0712	#0742	#0762
1/4"	#0713	#0743	#0763

DISS Female x DISS Female



Oxygen	Air	Vac
#0104	#SP-88	#SP-89

DISS Male x DISS Female



Oxygen	Air	Vac
#SP-95	#SP-96	#SP-97

DISS Female x 1/8" NPT Male



Oxygen	Air	Vac
#0110	#0209	#0309

DISS Male x 1/8" NPT Female



Oxygen	Air	Vac
#0702	#SP-98	#SP-99





50 PSI
Auxiliary
Valve

Part # 0150
(OXYGEN)

Part # 0151
(AIR)


Thread Comparisons Guide:



(Measurements given are approximate.)

For Oxequip, Medstar and Schrader
quick connects

PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM



NEB DRIVE RESCUE FLOW

An Economical adapter that delivers a preset flow from standard wall outlets.
It's a convenient, dependable way to administer aerosols.

DISS FEMALE



O.	6 to 7 LPM	8116
O.	7 to 8 LPM	8117
O.	8 to 9 LPM	8118
O.	10 to 11 LPM	8119
Air	6 to 7 LPM	8216
Air	7 to 8 LPM	8217
Air	8 to 9 LPM	8218
Air	10 to 11 LPM	8219
Rescue Flow	14 to 16 LPM	8115
O.	24 to 26 LPM	8114

OHMEDA



O.	6 to 7 LPM	8126
O.	7 to 8 LPM	8127
O.	8 to 9 LPM	8128
O.	10 to 11 LPM	8129
Air	6 to 7 LPM	8226
Air	7 to 8 LPM	8227
Air	8 to 9 LPM	8228
Air	10 to 11 LPM	8229
Rescue Flow	14 to 16 LPM	8125
O.	24 to 26 LPM	8124

CHEMETRON



O.	6 to 7 LPM	8136
O.	7 to 8 LPM	8137
O.	8 to 9 LPM	8138
O.	10 to 11 LPM	8139
Air	6 to 7 LPM	8236
Air	7 to 8 LPM	8237
Air	8 to 9 LPM	8238
Air	10 to 11 LPM	8239
Rescue Flow	14 to 16 LPM	8135
O.	24 to 26 LPM	8134

PURITAN-BENNETT



O.	6 to 7 LPM	8146
O.	7 to 8 LPM	8147
O.	8 to 9 LPM	8148
O.	10 to 11 LPM	8149
Air	6 to 7 LPM	8246
Air	7 to 8 LPM	8247
Air	8 to 9 LPM	8248
Air	10 to 11 LPM	8249
Rescue Flow	14 to 16 LPM	8145
O.	24 to 26 LPM	8144

**FOR CURRENT
PRICING INFORMATION
GO TO:
WWW.RCMEDICAL.COM**

REPLACEMENT RESPIRATORY OXYGEN SENSORS

Replaces: Ohmeda
6050-0004-110



Part # SS-20A

OEM Compatible with:
Analytical Ind.: PSR-11-915-4
C&S: G0-100
Ceramtec: CAG-10
City Technology: MOX-7
Envitec: OOM110
Datex-Ohmeda: 6050-0004-110
Integral Process: 65008
Maxtec: MAX-10, R112P10
Mercury Medical: 10-103-10
Nuova: E-15/4
Pacifitech: PT-20A
Sensidyne: SV-20A
Teledyne: R-30MED
Ventrex: V-20A
Vlamed: R-30V
VTI: 103105



Part # SS-03

OEM Compatible with:
Analytical Ind.: PSR-11-77
C&S: G0-070
Ceramtec: CAG-17
City Technology: MOX-4
Datex-Ohmeda: 6051-0000-219
Envitec: OOM105
Hamilton: HM-03
Hudson RCI: V-03, 5558
Integral Process: 65013
Maxtec: MAX-17, R116P10
Mercury Medical: 10-103-02
MSA: 804674
Newport: OM 100, VM 200
Nuova: E-77/O
Sensidyne: SV-03
SensorTech: ST-03
Teledyne: T-7
Ventrex: V-03
Vlamed: T-7
VTI: 103102

Fits ALL MiniOx!



Part # SS-11A

OEM Compatible with:
Air Shields: 6735142
Alphamed: R-17
Analytical Ind.: PSR-11-917-J
Smiths Medical > BCI: R-17
BioMed Devices: M-2, M-10
Bird Products: 66030
C&S: G0-030
Ceramtec: CAG-13
Datex-Ohmeda: 6051-0000-217
 Draeger: 6735142
Envitec: OOM103-1, OOM202-1
Hamilton: HM-11
Hudson RCI: 5567, V-11
Invivo: 9445
Lifecare: PLV-100, 37-079
Maxtec: MAX-13, R115P10
Mercury Medical: 10-103-05
MSA: 406931
Oxitron: 3310
Pacifitech: PT-11A
Medtronic > Puritan Bennett: 4-020933-00
Respironics: 37-079
Sensidyne: SV-11A
SensorTech: ST-11
Teledyne: R-17
Vlamed: R-41V
VTI: 100580



Part # SS-25A

OEM Compatible with:
Analytical Ind.: PSR-11-917-J1
C&S: G0-250
Maxtec: MAX-250ESF, R115P85
Mercury Medical: 10-103-08
Nuova: E-25/0
Pacifitech: PT-25A
Sensidyne: SV-25A
VTI: 103104



Part # SS-12A

OEM Compatible with:
Analytical Ind.: PSR-11-917M
C&S: G0-040
Ceramtec: CAG-12
City Technology: MOX-9
Datex-Ohmeda: 6051-0000-222
Envitec: OOM202
Hamilton Medical: 396008, 396009, HM-12
Hudson RCI: 5566, V-12
Integral Process: 65016
IMT Medical: PF300
Maxtec: MAX-12, R109P45
Mercury Medical: 10-103-06
MSA: 472062, 470644
Nuova: E-17/M
Pacifitech: PT-12A
Medtronic > Puritan Bennett: 4-072214-00
Siemens: 9004979-E347E, 6419332-E380E
Sensidyne: SV-12A
SensorTech: ST-01
Spacelabs: 90513
Taama: YR049700
Teledyne: R-24 MED, B-62321
Ventrex: V-12A
Vlamed: R-24Med
VTI Vascular Technology: 103100

ALL SENSORS

Unconditional 12 month Warranty!
Exceptional Quality and Performance!
Canned for increased shelf life and freshness!



FOR MORE SENSORS
PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM

MINIOX 200E OXYGEN ANALYZER



MiniOX 200E Oxygen Analyzer


Part # 600210

Comes complete with sensor, cable, tee and AA batteries

Replacement Sensor # 606940

The MiniOX 200E is best used anywhere an inline analyzer is needed, when alarms are not deemed important. It is especially useful in conjunction with older ventilators or oscillators. Designed with simplicity and portability in mind, this analyzer quickly adapts for use in tubing applications and is ideal for checking oxygen percent accuracy in line or on concentrations. Its patented ergonomic design fits comfortably in your hand and features one-touch on/off and calibration. The MiniOX 200E will calibrate to room air (20.9%) or 100% Oxygen. It is ideal for applications needing an external sensor.

- ① Large LCD – Easy-to-read characters
- ② Calibration Reminder – Reminds once a week
- ③ One-Touch Calibration – Auto calibrate to room air or 100% Oxygen
- ④ On/Off Button – One-touch on and off
- ⑤ Patented Ergonomic Design – Fits comfortably in your hand
- ⑥ External Sensor Cable – 10 ft. extended, 3 ft. retracted
- ⑦ Low Battery Indicator – 5,000+ Hours of battery life
- ⑧ Rubber Gasket – Prevents moisture from contact board
- ⑨ External Sensor – MiniOX 250E, 0-100% Oxygen, 24 month warranty



Approximately 5,000 hours with continuous use

MISCELLANEOUS

PLASTIC BAGS (DISPOSABLE)



PATIENT SET-UP (Cotton Drawstring)

Part No.	Description	Quantity
0049	8 x 10	500
0059	11 3/4 x 15	500
XL-S2	16 x 18 x 3	250

EQUIPMENT COVERS

Dust covers designed to protect ventilators and other equipment. Comes on a roll.

Part No.	Description	Quantity
PLC-60	28 x 22 x 58 (clear)	55
PLC-60/B	28 x 22 x 58 (buff)	55
0084	28 x 22 x 56 (clear***)	50
7083	28 x 22 x 56 (yellow)	50
0126	48 x 25 x 42 (red)	50

***printed, "CLEAN, READY FOR USE"

Main Flow Bacterial / Viral Filter



- 22mm ID x 15mm ID/22mm OD
- Bacterial filtration efficiency 99.999+%
- Viral filtration efficiency 99.99+%
 - Dead space: 42mL
- Flow resistance: approximately 1.5 cmH₂O at 60 LPM
- Individually packaged 50 to a case

Part Number 1605
Case of 50

ANESTHESIA ASSOCIATES PRESSURE MANOMETER



Part#	Description
199	Gauge with 1/4" NPT
274	Gauge with 5-55 cm H ₂ O adjustable pressure relief valve and "T" for supply tubing.



Part#	Description
199-MAX	Gauge with Pointer
274-MAX	274 with Pointer

Manometer has range of -40 to +80 cm H₂O in a 2 1/2" diameter chrome plated case with 1/4" NPT bottom connection. Has front zero correcting screw.

MAX is a resetable low friction pointer to capture maximum pressure or vacuum reading

POSEY CUFFLATOR™

ENDOTRACHEAL TUBE INFLATOR & MANOMETER



#8199 Cufflator with instrument case

An accurate, convenient tracheal cuff inflator and manometer. No more syringes, electricity or complicated connections needed to check cuff pressure regularly. Air vent button and inflator bulb let you vary cuff pressure quickly and easily. Clip on back allows easy positioning. Gauge shows recommended pressure range in cm H₂O.

**FOR MORE PRODUCTS AND SPECIALS
PLEASE VISIT OUR WEBSITE
WWW.RCMEDICAL.COM**

#2200 Flexible Connector (Silicone Rubber)

22mm I.D. x 22mm I.D.

This adapter was designed almost identical to the ones listed below:

INTEC 2200 MARQUEST 172200
BENNETT 3443 AM SYSTEMS 1139
Instrumentation Ind. 101-00-3443



Part Number 2200
Case of 50

Part Number 2200E
Eaches

MISCELLANEOUS

Timeter Airdyne Hospital Grade Air Compressors



Uniquely designed for bed-side use. Special material absorbs most noise for quiet performance. Audible and visual high temperature and low pressure alarms with DC powered six-volt backup. Vibration-free gauge shows diminished pressure, impaired drying or exceeded flow capacity. Circuit breaker and thermal overload protector included. Continuous flow of 64 LPM at 50 psi.

Part Number TE20



Designed for bedside use. Quiet, vibration-free performance. Provides condensation-free compressed air to respiratory therapy equipment. Vibration-free gauge shows diminished pressure, impaired drying or exceeded flow capacity. Circuit breaker and thermal overload protector included. Continuous flow of 42 LPM at 50 psi.

Part Number TE35

DeVilbiss Pulmo-Aide® Medication Nebulizer



3 PRONG HOSPITAL GRADE PLUG

Durable AC-powered medication nebulizer with Pulmo-Neb® 60-day disposable nebulizer and tubing. Unit weighs 7 lbs. and has 9 - 16 LPM range. 5 year warranty. **115 vac**

DV4650D

BLENDERS - PRECISION MEDICAL

Air/O₂



Its unique. It is the only unit available that allows the clinician to turn the bleed control to ON or OFF providing quieter operation. This feature is especially useful in the NICU and PICU.

Part #	Flow	Rate
5200	High Flow	15 - 120 lpm
5300	Low Flow	3 - 30 lpm

Heliox



Its unique. It is the only free-standing pneumatic HeliO₂ blending device on the market today. The ability to turn the bleed control collar ON or OFF provides ± 3% F_{IO2} accuracy as well as quieter operation and maximal Heliox tank durations.

Part #	Flow	%
5570	High Flow	70 - 30
5580	High Flow	80 - 20
5470	Low Flow	70 - 30
5480	Low Flow	80 - 20

PRECISION MEDICAL EASY AIR COMPRESSOR



EasyAir Compressor is a high performance portable medical air compressor designed to reliably supply compressed air 24 hours a day. EasyAir is ideal for continuous use in any clinical environment. Its sound dampened case is easy to clean, allows quiet operation and is durable enough for years of heavy use.

PM15P Pressure Gauge
PM15F Flow Gauge

**RC MEDICAL SELLS A WIDE VARIETY OF
RESPIRATORY SUPPLIES AND EQUIPMENT.
PLEASE VISIT OUR WEBSITE AT:
WWW.RCMEDICAL.COM**

REGULATOR AND FLOWMETER ACCESSORIES




Nut and Nipple (Christmas Tree)

Part #	Description	Quantity
7-4A*	One piece Green	Case/50
7-5A*	Two piece Green	Case/50
7-5C*	Two piece Yellow	Case/50
7-5H*	Two piece Clear	Case/50

Plastic - For single Patient Use (Above)*

0159	Metal Nut/Nipple	Each
0159C	Metal w/ 6" chain	Each

CYLINDER WRENCHES FOR SMALL AND LARGE CYLINDERS




Part #	Description	Quantity
7-92A	Plastic Small Wrench, BLACK	Case/10
7-92B	Plastic Small Wrench, GREEN	Case/10
7-92D	Plastic Small Wrench, PINK	Case/10
66071	Plastic Small Hand Wheel	Case/10
MCW-2B	Metal Small Cylinder Wrench	Case/10
MCW-2BC	Metal Small Wrench with Chain	Case/10
WR-3	Metal Wrench Small & Large	Each



REPLACEMENT GAUGES FOR OXYGEN REGULATORS

2" Bottom Port 1/4" NPT



1 1/2" Back Port 1/8" NPT



Low Profile 3000 PSI



Part #	Description
2510	15 lpm - 2" Bottom Port
2530	Contents (0-4000 PSI) 2" Bottom
2543	100 PSI 2" Bottom Port
2501	15 lpm - 1 1/2" Back Port
2521	Contents (0-4000 PSI) 1 1/2" Back
2500	15 lpm - 1 1/2" Bottom Port
2520	Contents (0-4000 PSI) 1 1/2" Bottom
MJ-C1	Low Profile Contents (3000 psi)
AIR, CO2 GAUGES	
MG-C11	AIR, CO2 (0 - 100 PSI) 2" Bottom
MG-C9	AIR, CO2 (0-4000 PSI) 2" Bottom
MG-C13	AIR, CO2 (0-4000 PSI) 1 1/2" Bottom

Yoke Washer

For Small Cylinder Regulators (CGA 870)



Part #	Description	Quantity
7-53E	aluminum w/ rubber center	each
7-53C	aluminum w/ rubber center	Case/50

LITERMETER

Every therapist should have one for checking the accuracy of regulators and flowmeters.



Part #	LPM
180-5004	0-2.5
180-5005	0-8
180-5006	2-15

ALL IN STOCK ORDERS PLACED BY 2 PM (EST) SHIPPED THE SAME DAY.

"Y" BLOCKS

Ohmeda Style

(Ohio)



2152	O ₂
2252	Air
2352	Vac

Chemetron Style

(NCG)



3152	O ₂
3252	Air
3352	Vac

REUSABLE NUT AND NIPPLE ADAPTER (CHRISTMAS TREE)

**HEX NUT
TUBING NIPPLE**



0159

0159
0159C

**HAND TIGHT
TUBING NIPPLE**



0154

0154
0154C

**60cm H₂O VACUUM
PRESSURE GAUGE**

- Reusable
- Read 60cm H₂O in 2cm Graduations
- Bottom 1/8" Small-bore Tubing Connection
- Packaged Individual



Part Number: NS60-PBS

**YOKE WASHER
(FOR SMALL CYLINDER
REGULATORS)**



7-53E Each
 7-53C Case/50

**ADJUSTABLE HEIGHT
SMALL CYLINDER CART
FOR D & E CYLINDERS**



#PX-1060